PaX-i™

vatech PaX-i™ User manual

PaX-i™

User Manual

Model : PCH-2500
Version : 2.88

- ENGLISH
- FRANÇAIS
- PORTUGUÉS
- ESPAÑOL


For the American regions

vatech

Postal code: 18449
13, Samsung 1-ro 2-gil, Hwaseong-si, Gyeonggi-do, Korea
www.vatech.com

vatech

ENGLISH

1. PaX-i Imaging System Overview.....	06
2. Getting started	08
3. Acquiring PANO Images (QR Code).....	13
4. Acquiring CEPH Images (QR Code).....	32

Notice

Thank you for purchasing the **PaX-i** extra-oral imaging system.

The **PaX-i** is an advanced digital dental diagnostic system that incorporates a Panoramic and Cephalometric (optional) imaging capabilities into a single system.

This manual describes how to operate the PaX-i system. It is recommended that you thoroughly familiarize yourself with this manual in order to make the most effective use of this equipment. Observe all cautions, safety messages and warnings which appear in this manual.

The **QR** code linked to video file about image capture for the **PANO Standard** and **CEPH Lateral** is provided in the manual. The smart phone and pad which have the QR code reader application program can be used to watch video.

Due to a constant technological improvement, the manual may not contain the most updated information, subjecting to change without prior notice to the persons concerned. For further information not covered in this manual, please contact us at:

Tel: +82-1588-9510

Email: gcs@vatech.co.kr

Website: www.vatech.com

Manual Name: PaX-i (PCH-2500) User Manual

Version: 2.88

Publication Date: 2023-06

Copyright © 2018 by VATECH

All rights reserved.

The documentation, brand name and logo used in this manual are copyrighted.

No part of this manual may be reproduced, transmitted, or transcribed without the expressed written permission of the manufacturer.

We reserve the right to make any alterations which may be required due to technical improvement.

Conventions in this Manual

The following symbols are used throughout this manual. Make sure that you fully understand each symbol and follow the instructions which accompany it.

	WARNING	indicates information that should be followed with the utmost care. Failure to comply with a warning may result in severe damage to the equipment or physical injury to the operator and/or patient.
	CAUTION	indicates a situation that demands prompt and careful action, a specific remedy, or emergency attention.
	IMPORTANT	indicates a situation or action that could potentially cause problems to the equipment and/or its operation.
	NOTE	emphasizes important information or provide useful tips and hints.
	X-RAY	indicates a possible danger of exposure to radiation.

1. PaX-i Imaging System Overview

The **PaX-i** is an advanced digital dental diagnostic system that incorporates panoramic and cephalometric (optional) imaging capabilities into a single system.

The **PaX-i** system is intended for use to take a panoramic and cephalometric images of the oral and craniofacial anatomy to provide diagnostic information for adult and pediatric patients.

The **PaX-i** system uses the advanced sensors to produce the higher quality of image in 2D of the head region, including the dental/maxillofacial regions, for planning and diagnostics.

The system includes the following software packages:

- **EasyDent** for 2D viewing and patient management
- **Imaging software**

The PaX-i may only be used by dentists, X-Ray technicians and other professionals who are licensed to perform X-Rays by the law of the region in which it is being used. (The PaX-i may be used for ENT (Ear, Nose and Throat) examination in countries requiring the CE marking.)

Standards and Regulations

The PaX-i is designed and manufactured to meet the following standards:
 MEDICAL - APPLIED ELECTROMAGNETIC RADIATION EQUIPMENT
 AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL
 HAZARDS ONLY IN ACCORDANCE WITH
 ANSI/AAMI ES 60601-1:2005 + AMD1:2012 + AMD2:2021
 CAN/CSA-C22.2 No. 60601-1:14 (Amendment 2:2022)
 IEC 60601-1-3:2008 + AMD1:2013 + AMD2:2021
 IEC 60601-2-63:2012 + AMD1:2017 + AMD2:2021
 NEMA standard publication PS 3.1-3.18, 2008

The CE symbol grants this equipment compliance with the European Directive for Medical Devices 93/42/EEC as amended by 2007/47/EC as a class IIb device.

MEDICAL - APPLIED ELECTROMAGNETIC RADIATION EQUIPMENT
 AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL
 HAZARDS ONLY IN ACCORDANCE WITH
 ANSI/AAMI ES 60601-1:2005 + AMD1:2012 + AMD2:2021
 CAN/CSA-C22.2 No. 60601-1:14 (Amendment 2:2022)
 IEC 60601-1-3:2008 + AMD1:2013 + AMD2:2021
 IEC 60601-2-63:2012 + AMD1:2017 + AMD2:2021


Classifications (IEC60601-1 6.1):

Protection against the ingress of water: Ordinary Equipment: IPX0

Protection against electric shock: Class 1 equipment, Type B Applied Parts (Chinrest, Bite and cover, Nasal Positioner and cover, Ear Rod and cap, Carpus Plate)

General View of the PaX-i

NO.	Item	NO.	Item
1	X-RAY SENSOR (CEPH)	11	CHIN SUPPORT
2	EAR ROD	12	COLUMN UP / DOWN BUTTON 
3	NASAL POSITIONER	13	LASER BEAM ON / OFF BUTTON 
4	SECONDARY COLLIMATOR	14	EMERGENCY SWITCH
5	BASE	15	HANDLE FRAME
6	LED LAMP	16	ON / OFF SWITCH
7	ROTATING UNIT	17	TELESCOPIC COLUMN
8	X-RAY SENSOR (PANO)	18	COLUMN UP / DOWN SWITCH (OPTIONAL)
9	X-RAY TUBE	19	STATIONARY COLUMN
10	TEMPLE SUPPORT		

2. Getting started

2.1 Turning on the PaX-i


Do not position the patient in the unit while it is initiating. The patient could be injured if the equipment malfunctions.

- Condensation could form inside the equipment if it is a different temperature than the surrounding room. Only turn on the equipment once it has reached room temperature.
- Wait at least 20 seconds after the equipment has been turned off to turn it back on.
- Allow the equipment to warm up for at least 5 minutes before acquiring images or preferably more than 30 minutes for image quality.

Turn the **ON / OFF** switch underneath the handle frame to the ON position. Check whether the LED lamp on top of the equipment is illuminated. If it is illuminated, the equipment is ready for image acquisition.


2.2 Running the Image Viewer – EasyDent

EasyDent is a basic imaging platform for all VATECH's dental X-ray equipments. Both the **Imaging Program** and **3D Viewer** are interfaced with **EasyDent**.

On your desktop, double-click **EasyDent** or click **Start** → **All Programs** → **EasyDent**.

The **EasyDent** main window will be displayed.

Go to the **section 2.2.1 ~2.2.2** to create or open an existing patient record.


NOTE

For detailed guide to using this software, please refer to the **EasyDent user manual**.

2.2.1 Creating a New Patient Record

To create a new patient record, follow the procedure outline below:

1. Click the **Patient icon** () on the upper left corner of the EasyDent's main GUI window.

2. The following dialog box will open.

The 'Add Patient' dialog box contains the following fields and controls:

- Chart No.: 20111223_180547 (with 'Auto' and 'Auto No.' checkboxes)
- First Name: [Empty text box]
- Last Name: [Empty text box]
- Social ID: [Empty text box]
- Birthday: 2011 (Year), 1 (Month), 1 (Day)
- Gender: Male (Dropdown)
- Treatment: Treatment (Dropdown)
- Address1: [Empty text box]
- Address2: [Empty text box]
- E-mail: [Empty text box]
- Tel: [Empty text box]
- Mobile: [Empty text box]
- Doctor: [Empty dropdown menu]
- Buttons: Add, Cancel

3. Enter the required patient information. **Chart Number, First Name, and Last Name** are required fields which must be filled in. All other fields are optional, but it is recommended that they be filled in.

4. Click **Add** to save the patient record.


2.2.2 Retrieving Patient Records

You can search through the patient database using a patient's chart number, their first name, or their last name.

1. On the **Patient information pane**, double-click the **Chart No., First name, or Last name** of the patient and the virtual keyboard will pop up.


2. Enter the **Chart No.**, **First name**, or **Last name** of the patient by clicking the mouse on the virtual key board and click the **Enter**. (The physical keyboard can be used to do the same job.)
3. Patient information can be displayed on the Patient information pane and Patient List.

2.3 Initiating the Imaging Program

1. First, click the patient information in the patient list, and click the **PANO/CEPH** icon in the upper left corner of the EasyDent's main window to open the imaging program.

2. The following imaging program window opens. The sole purpose of this window is to control equipment settings and acquire images.

Please proceed to the next chapter.


NOTE

Refer to chapters 3 ~ 4 for information regarding image acquisition.

3. Acquiring PANO Images

3.1 Setting Exposure Parameters

Perform the following procedures to select the capture parameters for the specific patient and capture mode.

1. Choose an imaging mode under **Imaging Mode**.

2. Select the mode under **PANO Type**.

Mode	Details	
Normal	Normal image	
Magic PAN	Image with ultra-high resolution	(optional)

3. Select the image quality for the image under **Image Quality**.


Mode	Details	
UHD	Image with Ultra-high resolution	(optional)
HD	Image with higher resolution than Normal mode Takes longer scan time than the Normal .	
Normal	Normal image	

4. Select the patient's arch type under Arch Selection. By default, Normal is selected.

Orthogonal Mode: This mode enables overlapping regions of teeth to be minimized when acquiring images in the ROI. If orthogonal arch is selected, its sub-modes are activated.

5. Select the ROI for panoramic image acquisition under **PANO Examination**.

6. The patient's gender will be automatically selected according to the patient information registered in EasyDent. Nevertheless, ensure that this information is correct.

Age Group		VATECH's Standard
	Child	≤ 12
Adult	Man	≥ 13
	Woman	

7. Select X-ray intensity.

X-ray intensity (Hard, Normal, Soft) is according to the operator's decision.

Soft ≤ Normal ≤ Hard

8. A default value for the tube voltage (kVp) and current (mA) will be displayed based on the gender and X-ray intensity. If necessary, further adjustments can be made using the arrows to the right of each number.

9. Click the **Confirm** button for these parameters to take effect.

Please wait for a moment, while the rotating moves to its initial scanning position.

When you click the **Confirm** button,

- The **Ready** button will begin blinking to show that it has been activated. This means that the equipment is ready for X-ray exposure.
- The rotating unit will move to its initial scanning position.
- Three laser beams, **Mid-sagittal plane, Frankfurt plane and Canine**, will be activated to facilitate patient positioning. These beams will disappear after 20 minutes or if the **Ready** button is clicked.

The laser beam on/off icons are located:

- In the upper right side of the window: 
- On the Handle Frame: 

Scan Time : 0.0
DAP : 0.00 uGy x m²


- The scan time and estimated DAP (Dose Area Product) values will be shown on the main display for the selected exposure.

10. Guide the patient to the equipment. Position the patient within the equipment. For further information about patient positioning, refer to section **3.2. Positioning the Patient**.


3.2 Positioning the Patient

Before Positioning the Patient

- Ask the patient to remove all jewelry and metallic objects, such as earrings, hair pins, spectacles, dentures, and orthodontic appliances. These items may cause shadow images which could obscure diagnosis.
- It is strongly recommended that the patient wear a lead apron for protection against any possible scatter radiation.
- Adjust the column height to the patient's height using the up/down button or switch (optional).

Column up/down buttons on the handle frame

Column up/down switch (optional)

NOTE

Correct positioning is an important factor in capturing the best possible image. Proper positioning will reduce the appearance of the cervical spine in the image.

CAUTION

If you don't adjust the beams for patient alignment (vertical, horizontal, and canine beams), part of the image may be zoomed in or out or a ghost may appear in the image.

WARNING

Be careful not to project the laser beams directly into the patient's eyes as this could severely damage the patient's vision.


3.2.1 PANO Standard Mode

Positioning the patient


PANO Standard_ENG

*To link to the video manual: scan the QR code with smart phone or pad.

1. Insert the normal chin support and bite block into the chin support receptacle on the equipment.

2. Place a hygienic cover over the bite block.

*Chin support: Edentulous

WARNING


The hygiene cover for the bite block is intended for single use only. Always replace the hygiene cover for each new patient.

CAUTION

Disinfect the chin support and bite block using an alcohol-based cleaning solution and wipe away all residue with dry cloth.


Temple support wheel


3. Loosen the Temple Supports by turning the Temple Support Wheel.

4. Guide the patient — facing the chin support — to the equipment.

Column up / down
button


5. Adjust the height of the column using the column up/down button or switch (optional) until the patient's chin rests comfortably on the chin support.

6. Ask the patient to:

- Stand upright
- Firmly grip the handle
- Lean his/her chest lightly against the equipment
- Position his/her feet slightly forward and open stance

7. Ensure that the patient's shoulders remain level and that his/her neck is relaxed. The cervical spine should be straight and upright.

8. Have the patient bite the bite block along its grooves using his/her upper and lower incisors (for edentulous patients, use the edentulous chin support).

9. Have the patient

- Close his/her lips around the bite block
- Keep his/her tongue pressed against his/her palate
- Close his/her eyes

Ask the patient to remain still until scanning is completed.

CAUTION

To acquire the best image possible, ask the patient not to:

- Breathe or swallow saliva during image acquisition
- Move during image acquisition

Aligning the Laser Beams


WARNING

Be careful not to project the laser beams directly into the patient's eyes as this could severely damage the patient's vision.


1. **Mid-sagittal plane laser beam:** Position the Mid-sagittal plane laser beam on the center of the patient's face to prevent magnifications of the left or right side in the final image.

2. **Frankfurt plane laser beam:** Position the patient's head so that the Frankfurt plane is aligned with the Frankfurt plane laser beam.

** The Frankfurt plane is the plane which joins the infra-orbital point to the superior border of the external auditory meatus.*


3. **Canine laser beam:** Have the patient smile to properly position the canine laser beam on the center of the patient's canine tooth.

Completing Patient Positioning

1. After confirming patient positioning and beam alignment, adjust the temple supports to fit snugly on both sides of the patient's head using the temple support wheel. The temple support wheel is located underneath the patient support rest.

2. Click the **Ready** button on the GUI after the patient has been properly positioned. No X-Ray will be emitted at this point. The rotating unit will then move into position for imaging.

3. Proceed to section **3.3. Initiating X-Ray Exposure**.

3.2.2 TMJ Open Mode


There are two TMJ sub modes: TMJ Open and TMJ Close. Scan a TMJ Open image first and then scan a TMJ Close image. To prepare and position the patient, follow the process outlined below.

- **Correct positioning is an important factor in capturing the best possible image. Proper positioning will reduce the appearance of the cervical spine in the image.**

- **For a child with a large head circumference, *Man/Woman* mode may be selected instead of *Child*.**


Positioning the patient

1. Insert the TMJ chin support into the receptacle on the equipment.

Disinfect the chin support using an alcohol-based cleaning solution and wipe away all residues with dry cloth before proceeding any further.


Temple support wheel

2. Loosen the **temple Supports** by turning the **temple support wheel**.

3. Guide the patient — facing the chin support — to the equipment.
4. Adjust the height of the equipment using the column up/down button or switch (optional) until the patient's chin is resting on the chin support. Ensure that the chin is in contact with the chin support.
5. Ask the patient to:
 - Stand upright
 - Firmly grip the handle
 - Lean his/her chest lightly against the equipment
 - Position his/her feet slightly forward.

Patient position with chin support

6. For **TMJ Open** imaging, have the patient:

- Open his/her mouth as wide as possible
- Keep his/her tongue positioned towards the palate
- Close his/her eyes.
- Be careful the chin not to touch any part of the equipment.

The top of the TMJ support and the patient's acanthion point must remain in contact at all times during exposure.

7. Ensure that the patient's shoulders remain level and that his/her neck is relaxed. The cervical spine should be straight and upright.

8. Ask the patient to remain still until scanning is complete.

To acquire the best image possible, ask the patient not to:


- Breathe or swallow saliva during image acquisition
- Move during image acquisition

Aligning the Laser Beams



Be careful not to project the laser beams directly into the patient's eyes as this could severely damage the patient's vision.

1. **Mid-sagittal plane laser beam:** Position the Mid-sagittal plane laser beam at the center of the patient's face to prevent magnifications of the left or right side in the final image.


2. **Frankfurt plane laser beam:** Position the patient's head so that the Frankfurt plane is aligned with the Frankfurt plane laser beam.

3. **Canine laser beam:** Have the patient smile to properly position the canine laser beam on the left corner of the patient's lips.

Completing Patient Positioning

1. After confirming patient positioning and beam alignment, adjust the temple supports to fit snugly on both sides of the patient's head using the temple support wheel. The temple support wheel is located underneath the patient support rest.

2. Click the **Ready** button on the GUI after the patient has been properly positioned. No X-Ray will be emitted at this point. The rotating unit will then move into position for imaging.


3. Proceed to section **3.3. Initiating X-Ray Exposure**.

3.2.3 TMJ Close Mode

Scan a **TMJ Open** image first and then scan a **TMJ Close** image. Guide the patient to close his/her lips for **TMJ Close** imaging while maintaining the same position used for **TMJ Open** imaging. The positions are the same for both modes, except for the position of the lips.

1. Shortly after **TMJ Open** imaging has been completed, you will receive the following message: "Do you want to take the **TMJ Close** image continuously?" Press the **OK** button to capture a **TMJ Close** image.

2. For **TMJ Close** imaging, Have the patient:

- Close his/her mouth
- Keep his/her tongue pressed against his/her palate
- Close his/her eyes


The top of the **TMJ** support must remain in contact with the patient's acanthion point during exposure.

Ask the patient to remain still until scanning is complete.

3. The laser beams should be aligned in the same way as for **TMJ Open** imaging.
4. Proceed to section **3.3. Initiating X-Ray Exposure**.

3.2.4 Sinus Mode


Positioning the Patient

1. Insert the chin support (Sinus) into the receptacle on the equipment.


Disinfect the chin support using an alcohol-based cleaning solution and wipe away all residues with a dry cloth before proceeding any further.

2. Loosen the temple supports by turning the temple support wheel.


3. Guide the patient —facing the chin support — to the equipment. Have the patient stand upright in the center of the equipment.

4. Adjust the height of the equipment using the column up/down button or switch (optional) until the patient's chin rests comfortably on the chin support.


5. Ask the patient to:

- Stand upright
- Firmly grip the handle
- Lean his/her chest lightly against the equipment
- Position his/her feet slightly forward.

6. Ask the patient to press his/her lips against the chin support.

7. Ensure that the patient's shoulders remain level and that his/her neck is relaxed. The cervical spine should be straight and upright.

8. Have the patient:

- Tilt his/her head back by 10-15°
- Close his/her mouth
- Keep his/her tongue pressed against his/her palate
- Close his/her eyes

Ask the patient to remain still until scanning is completed.


To acquire the best image possible, ask the patient not to:

- Breathe or swallow saliva during image acquisition
- Move during image acquisition

Aligning the Laser Beams


Be careful not to project the laser beams directly into the patient's eyes as this could severely damage the patient's vision.

1. **Mid-sagittal plane laser beam:** Position the Mid-sagittal plane laser beam at the center of the patient's face to prevent magnifications of the left or right side in the final image.


2. **Frankfurt plane laser beam:** The Frankfurt plane laser beam should go from the top of the ear to the tip of the nose when the patient's head is tilted back 10 – 15 °.

3. **Canine laser beam:** Have the patient smile to properly position the canine laser beam on the center of the patient's canine tooth.

Completing Patient Positioning

1. After confirming patient positioning and beam alignment, adjust the temple supports to fit snugly on both sides of the patient's head using the temple support wheel. The temple support wheel is located underneath the patient support rest.

2. Click the **Ready** button on the PC after the patient has been properly positioned. At this point, X-Ray will not be emitted. The rotating unit will move to position for imaging.

3. Proceed to section **3.3. Initiating X-Ray Exposure**.


3.3 Initiating X-Ray Exposure

The X-Ray exposure method and its sequences are the same for PANO Standard and PANO Special modes. The example used in this manual is the PANO Standard mode.

To emit the X-Ray, perform the following procedure.

If a problem occurs during image acquisition, press the red emergency stop button to immediately stop all moving parts and cut off all power to the equipment's electrical components. You may then safely release the patient from the equipment.

To reset this button, turn it clockwise until it pops up.


- Do not operate the PC during exposure. Failure to comply with this instruction may cause the system to malfunction.
- The operator shall observe the X-ray safety regulations applicable to their area at all times during the operation of this equipment.

If an emergency occurs during image acquisition, release the exposure switch to cease X-ray emission.

1. Leave the X-Ray room and close the door. The operator must visually monitor the patient at all times during image acquisition.
2. Press and hold down the exposure switch until image acquisition is complete.


During X-Ray irradiation, ensure that

- The LED lamp on top of the equipment turns orange to indicate X-Ray emission.
- The warning lamp outside the X-Ray room turns on.
- The sound (beep or music: optional) goes off.
- The radiation symbol on the upper left corner of GUI turns yellow to indicate X-Rays are being emitted.

3. The image appears in real time on the imaging GUI.

CASE 1: Normal / AF

Normal / AF: Ex) 5 images view is set

To display an image at a time on the screen, either one can be used: clicking the radio button or scrolling the mouse wheel. When clicked on the **Save** button, the current image is saved.


NOTE

With AF feature enabled, the best image of the 2 or 5 images taken can be saved. This option can be selected by the user from the settings screen by following the process below: Settings () → PANO / CEPH → Multi focusing setting: Select 2 or 5 images → Save.


If, for a single image view, the automatic save as a default is selected, the captured image will be saved automatically.

CASE 2: Magic PAN


Two images are obtained, with each in the **Magic PAN** and **Normal** respectively. You would click either button to compare them with each other.

The buttons:

4. Click the **Save** button to save the image.

5. The captured image will be automatically transferred to **EasyDent**. Clicking the patient's name in the patient list after imaging will refresh the image list for that patient. The patient's most recent image will appear on the far left of the thumbnail image pane, as shown below.

6. Double-click an image to enlarge it for viewing or to check the quality of the image.

After Image Acquisition


After acquiring the image, complete the following steps:

- Loosen the temple supports to release the patient.
- Remove the hygiene barrier from the bite block (for standard PANO mode only).


4. Acquiring CEPH Images

4.1 Setting the Exposure Parameters

Perform the following procedures to set the exposure parameters for the specific patient and capture mode. Depending on the sensor type employed, one of two kinds of imaging S/W comes with the equipment for the CEPH mode examination.

With one shot type sensor

With scan type sensor

1. Click the **CEPH** button.
2. Select the scanning mode under **CEPH Examination**. – **Lateral, PA, SMV, Waters View, Carpus**.
3. **Exposure time:** If necessary, this time can be adjusted by resolution of 0.1 s in the range of 0.7 s to 1.2 s.

4. The gender and age of the patient is automatically selected according to EasyDent's patient information. Nevertheless, ensure that this information is correct.

Age Group		VATECH's Standard
	Child	≤ 12
Adult	Man	≥ 13
	Woman	

NOTE

A Child is defined as a person who is younger than 12 years old. If child is selected, the image size and exposure dose are automatically reduced.

5. Select X-ray intensity.

NOTE


X-ray intensity (Hard, Normal, Soft) is according to the operator's decision.

Soft ≤ Normal ≤ Hard

6. A default value for the tube voltage (kVp) and current (mA) will be displayed based on the gender and X-ray intensity.

If necessary, you may manually adjust these values using the arrows to the right of each value.

7. Click **Confirm** to accept the parameters.

When you click the **Confirm** button,

- The **Ready** button will begin blinking to show that it is activate. This means that the equipment is ready for X-Ray exposure.

- The CEPH sensor and the 2nd collimator enter positioning mode.

Scan Time : 0.0
DAP : 0.00 uGy x m²

- The scan time and estimated **DAP** (Dose Area Product) value are shown on the main display for the intended exposure.


8. Guide the patient to the equipment. Position the patient. Refer to section **4.2. Positioning the Patient**.

4.2 Positioning the Patient


Follow the procedure below to prepare and position the patient for CEPH imaging.

Before Positioning the Patient

- Ask the patient to remove all jewelry and metallic objects, such as earrings, hair pins, spectacles, dentures, and orthodontic appliances. These items may cause shadow images that could impede diagnosis.
- It is strongly recommended that the patient wear a lead apron for protection against any possible scatter radiation.
- Adjust the column height to the patient's height using the up/down button or switch (optional).

Column up/down buttons on the handle frame

Column up/down switch (optional)

NOTE

Correct positioning is an important factor in capturing the best possible image.

Ensure that the nasal positioner left unfolded, before adjusting the ear rods in the proper direction.

Ensure that the patient is clear of all moving parts while adjusting the height of the equipment.


4.2.1 Lateral Mode

CEPH Lateral_ENG

*To link to the video manual: scan the QR code with smart phone or pad.

Follow the procedure below to position the patient.


1. Widen the distance between the two ear rods.

Porion position reference indicator enables the operator to easily identify the position of the porion on the image.

2. Guide the patient to the CEPH unit.


3. Adjust the height of the unit to suit the patient by pressing the column up/down button or switch (optional).

4. Ask the patient to stand upright.

5. Make sure that the patient's shoulders are level and that their neck is relaxed.

6. The ear rods should fit properly along the patient's ear canals. The patient's Frankfurt Plane should be parallel to the floor.


7. Position the nasal positioner on the patient's Nasion point. The height of nasal positioner can be adjusted.

After adjusting the height of the column to suit the patient, fit the ear rods along the patient's ear canals and adjust the nasal positioner.

8. Ask the patient to swallow saliva and remain still with his/her mouth closed until X-Ray exposure is complete.


9. Click the **Ready** button after the patient has been properly positioned. No X-Ray will be emitted at this point. The CEPH sensor and the 2nd collimator will enter positioning mode.

10. Proceed to section **4.3. Initiating X-Ray Exposure.**

4.2.2 Frontal Mode (PA)

Follow the procedure below to correctly position the patient.


1. Turn the ear rods 90° clockwise from their initial position.

2. Widen the distance between the two ear rods.


NOTE

Porion position reference indicator enables the operator to easily identify the position of the porion on the image.

3. The nasal positioner should be flipped to the side and upward to prevent it from obstructing the image acquisition.

4. Guide the patient to the CEPH unit.

5. Adjust the height of the unit to suit the patient by pressing the column up/down button or switch (optional).


6. Ask the patient to stand upright facing the sensor.

Make sure that the patient's shoulders are level and that his/her neck is relaxed.


CAUTION

After adjusting the column to the patient's height, fit the ear rods along the patient's ear canals.

7. The ear rods should fit properly along the external auditory meatus of the patient. The patient's Frankfurt Plane should be parallel to the floor.

8. Ask the patient to swallow saliva and remain still with his/her mouth closed until X-Ray exposure is complete.


9. Click the **Ready** button after the patient has been positioned. No X-Ray will be emitted at this point. The CEPH sensor and the 2nd collimator will enter capturing mode.

10. Proceed to section **4.3. Initiating X-Ray Exposure**.


4.2.3 SMV Mode

Follow the procedure below to prepare and position the patient.

1. Widen the distance between the two ear rods.

2. The nasal positioner should be flipped to the side and up to prevent it from obstructing the image acquisition.

3. Guide the patient to the CEPH unit.


4. Adjust the height of the unit to suit the patient by pressing the column up/down button or switch (optional).

5. Ask the patient to stand upright facing the secondary collimator.

Make sure that the patient's shoulders are level and that his/her neck is relaxed.


After adjusting the column to suit the height of the patient, fit the ear rods along the patient's ear canals.

6. Place the ear rods along the patient's ear canals. Make sure that the ear rods are comfortably, yet firmly in place.

7. Gently tilt the patient's head back until his/her Frankfurt plane is perpendicular to the floor, as shown below.


8. Click the **Ready** button after the patient has been properly positioned. No X-Ray will be emitted at this point. The CEPH sensor and the 2nd collimator will enter positioning mode.

9. Proceed to section **4.3. Initiating X-Ray Exposure**.

4.2.4 Waters View Mode

Follow the procedure below to prepare and position the patient for **Waters View** mode image acquisition.

1. Widen the distance between the two ear rods.

2. The nasal positioner should be flipped to the side and up to prevent it from obstructing the image acquisition.

3. Guide the patient to the CEPH unit.

4. Ask the patient to stand upright facing the sensor.


5. Adjust the height of the unit to suit the patient by pressing the column up/down button or switch (optional).

6. Make sure that the patient's shoulders are level and that his/her neck is relaxed.


After adjusting the column to suit the height of the patient, fit the ear rods along the patient's ear canals.

7. Place the ear rods along the patient's ears. Make sure that the rods are comfortably, yet firmly, in place.

8. Ask the patient to swallow saliva and tilt his/her neck back 30° - 40°, with the mouth closed, until X-Ray exposure is complete.


9. Click the **Ready** button after the patient has been properly positioned. No X-Ray will be emitted at this point. The CEPH sensor and the 2nd collimator will move to capturing mode.


10. Proceed to section **4.3. Initiating X-Ray Exposure**.

4.2.5 CARPUS Mode

Attaching the CARPUS plate

1. Match the notch on the **CARPUS** plate to the end of nasal positioner.


2. Slide the **CARPUS** plate onto the nasal positioner.

3. Slide the locking bolt on the **CARPUS** plate towards the nasal positioner and turn it to secure the **CARPUS** plate in place.

4. Ensure that the **CARPUS** plate is firmly in place.

Positioning the Patient

1. Ask the patient to put his/her right hand flat on the **CARPUS** plate. It is important to ensure that the patient does not bend his/her fingers.

2. Have the patient close their eyes and remain still until scanning is complete.

Ask the patient not to block the ruler on the nasal positioner by his / her fingers. It could deteriorate image quality.

3. Click the **Ready** button after the patient has been properly positioned. No X-Ray will be emitted at this point.


4. Proceed to section **4.3. Initiating X-Ray exposure**.

4.3 Initiating X-Ray exposure

The X-Ray exposure method and process is the same for all CEPH modes. The example and images below are taken from an X-Ray performed in **CEPH Lateral** mode.

If a problem occurs during image acquisition, press the red emergency stop button to immediately stop all moving parts and cut off all power to the equipment's electrical components. You may then safely release the patient from the equipment.

To reset this button, turn it clockwise until it pops up.

Do not operate the PC during exposure. Failure to comply with this instruction may cause the system to malfunction.

The operator shall observe the X-Ray safety regulations applicable to their area at all times during the operation of this equipment.

If an emergency occurs during image acquisition, release the exposure switch to cease X-ray emission.


1. Leave the X-Ray room and close the door. The operator must maintain visual contact with the patient at all times during image acquisition.

2. Press and hold down exposure switch until the image is acquired.

Exposure indicator light turns orange

Orange : Exposure

X-Ray ON indicator


During X-Ray irradiation, ensure that

- The LED lamp on top of the equipment turns orange to indicate X-Ray emission.
- The warning lamp outside the X-Ray room turns on.
- The sound (beep or music: optional) goes off.
- The radiation symbol on the upper left corner of GUI turns yellow to indicate X-Rays are being emitted.

The image appears in real time on the imaging GUI.

3. Click the **Save** button to save the acquired image. If automatic save has been selected as a default setting, the image will be saved automatically.

4. The captured image will be automatically transferred to **EasyDent**. Clicking the patient's name in the patient list after imaging will refresh the image list for that patient. The patient's most recent image will appear on the far left of the thumbnail image pane.

5. Double-click on the image to enlarge it for viewing or to check its quality.


Lateral Mode with one shot type sensor
- Irradiated condition: 95 Kvp / 10 mA / 0.9 s

After Image Acquisition

After the image has been acquired, perform the following tasks:

- Fold away the nasal positioner.
- Loosen the ear rod supports and remove them from the patient's ears.
- Release the patient.

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Copyright © 2018 par VATECH
Copyright © 2018 por VATECH