vatech A9

vatech A9 User manual for the American region

vatech A9

User manual

Model : PHT-30CSS
Version : 1.10

- English
- Français
- Português
- Español

Notice

Thank you for purchasing the **vatech A9 (PHT-30CSS)** extra-oral imaging system.

vatech A9 (PHT-30CSS) is one of the manufacturer's product series for aiding dental professionals in providing excellent care in a safe environment that promotes healing.

vatech A9 (PHT-30CSS) is an advanced digital diagnostic system that incorporates PANO, CEPH (Optional), CBCT imaging capabilities into a single system.

This manual describes how to operate the **vatech A9 (PHT-30CSS)** system. It is recommended that you thoroughly familiarize yourself with this manual to make the most effective use of this equipment.

Observe all cautions, safety messages, and warnings in this manual.

Due to constant technological improvement, the manual may not contain the most updated information and is subject to change without prior notice to the persons concerned. For further information not covered in this manual, please contact us at:

VATECH Co., Ltd.

Phone: +82-1588-9510

E-mail: gcs@vatech.co.kr

This document is originally written in English.

vatech A9 (PHT-30CSS) is referred to as "equipment" in this manual.

Manual Name: vatech A9 (PHT-30CSS) User Manual

Version: 1.10

Publication Date: 2023-05

Document Number: VDH-UM-098

This Page is Left Blank Intentionally

Table of Contents

Notice.....	i
Table of Contents	iii
Conventions in this Manual	1
1. Imaging System Overview	3
1.1 System Components.....	3
1.2 Features.....	3
1.3 Imaging System Options.....	3
1.4 Standards and Regulations.....	4
1.5 Operating Principles.....	4
1.6 Equipment Overview	5
2. Getting Started.....	11
2.1 Turning on the Equipment.....	11
2.2 Running the Image Viewer (EzDent-i).....	12
2.3 Initiating the Console Software	14
3. Acquiring PANO Images	15
3.1 Configuring Exposure Parameters	15
3.2 Patient Positioning	19
3.3 X-ray Exposure	35
3.4 Finishing the Scan	36
3.5 Checking the Captured Images.....	36
4. Acquiring CEPH Images (Optional).....	37
4.1 Configuring Exposure Parameters	37
4.2 Patient Positioning	40
4.3 X-ray Exposure	52
4.4 Finishing the Scan	53
4.5 Checking the Captured Images.....	53
5. Acquiring CBCT Images.....	55
5.1 Configuring Exposure Parameters	55
5.2 Patient Positioning	59

Table of Contents

5.3	X-ray Exposure	65
5.4	Finishing the Scan	66
5.5	Checking the Captured Images	66
5.6	Image Comparison according to Patient Position – Incorrect Posture	67

Conventions in this Manual

The following symbols are used throughout this manual. Make sure that you fully understand each symbol and follow the instructions accompanied.

To prevent physical injury and damage to the equipment, please observe all warnings and safety information included in this document.

	WARNING	Indicates information that should be followed with the utmost care. Failure to comply with a warning may result in severe damage to the equipment or physical injury to the operator and patient.
	CAUTION	Indicates a situation that demands prompt and careful action, a specific remedy, or emergency attention.
	IMPORTANT	Indicates a situation or action that could potentially cause problems to the equipment and its operation.
	NOTICE	Emphasizes essential information or provides useful tips and hints.
	RADIATION	Indicates a possible danger from exposure to radiation.
	SINGLE USE	Indicates a component that must be replaced for each new patient.
	ESD susceptibility	Indicates that an item is susceptible to damage from electrostatic discharges.

This Page is Left Blank Intentionally

1. Imaging System Overview

1.1 System Components

- **vatech A9 (PHT-30CSS)** X-ray equipment
- PC system
- Console Software: PANO, CEPH (Optional), and CBCT Scan
- **EzDent-i**: 2D viewer and patient management software
- **Ez3D-i**: 3D viewer software

1.2 Features

- FOV 8x8 support (Anatomically 9.3x8.0) (cm)
- The multi-imaging solution for Accurate Diagnostics
- Conventional 2D (PANO and CEPH) image acquisition
- Control Panel implemented for easy use
- DICOM (Digital Imaging Communication in Medicine) format supported
- Differentiated Console Software Interface

1.3 Imaging System Options

Configuration	Item	Sensor	
SP	PANO +CBCT	PANO/CBCT	Xmaru1404CF-Plus
RC	PANO +CBCT +CEPH	PANO/CBCT	Xmaru1404CF-Plus
		CEPH	Xmaru2602CF

1.4 Standards and Regulations

Standards

vatech A9 (PHT-30CSS) is designed and developed to comply with the following international standards and regulations:

- MEDICAL - APPLIED ELECTROMAGNETIC RADIATION EQUIPMENT
ANSI/AAMI ES60601-1 (2005) + AMD 1 (2012)
CAN/CSA-C22.2 No. 60601-1 (2014)
IEC 60601-1-3:2008 + AMD1:2013, IEC 60601-2-63:2012 + AMD1:2017 21 CFR 1020.30, 31, 33
- NEMA Standard publication PS 3.1-3.18

	<p>This is Class IIb equipment and obtained CE marking in April 2007 for regulations compliance by the revised European Union's MDD (Medical Devices Directive) 93/42 EEC.</p>
	<p>MEDICAL - APPLIED ELECTROMAGNETIC RADIATION EQUIPMENT ANSI/AAMI ES60601-1 (2005) + AMD 1 (2012) CAN/CSA-C22.2 No. 60601-1 (2014) IEC 60601-1-3:2008 + AMD1:2013, IEC 60601-2-63:2012 + AMD1:2017</p>

Classifications (IEC 60601-1 6.1)

- The degree of protection against water ingress: Ordinary Equipment: IPX0
- The degree of protection against electric shock: Class 1 equipment, Type B Applied Parts (Chinrest, Bite and cover, Nasal Positioner and cover, Ear Rod and cap, Carpus Plate)


1.5 Operating Principles

X-ray is emitted when a high voltage is supplied to X-ray tube assembly, which frees electrons from the cathode.

They hit anode to produce an X-ray. The machine acquires images by emitting X-ray continuously and rotates on the human tooth at different angles.

Images are acquired, computed, and recompiled to reproduce 2D or 3D images.

1.6 Equipment Overview

<Perspective View>

No.	Item	Description
1	X-ray Detector for CEPH (Optional)	Xmaru2602CF for CEPH imaging sensor
2	Nasal Positioner	Positions the patient during CEPH imaging. The ruler used for reference in an acquired image that is different from the actual size
3	Ear Rods	Secure the patient's head during CEPH imaging.
4	Column UP/DOWN Switch (optional)	Adjust column height to adjust the height of the vertical frame.
5	Stationary Column	Supports the whole part of the equipment.
6	Base (Optional)	Balances the equipment and maintains its safety.
7	Main Power Switch	Turns on/off the main power of the equipment.
8	D-Sub Connector	The input signal port for Column UP/DOWN Switch
9	Emergency Stop Switch	Immediately stops the moving parts and cuts off all power to the equipment's electrical components.
10	Column UP/DOWN button	Adjust column height to adjust the height of the vertical frame.

1. Imaging System Overview

No.	Item	Description
11	Chinrest	The place to rest the chin.
12	X-ray Detector for PANO/CBCT	Xmaru1404CF-Plus for PANO/CBCT imaging sensor
13	LED Lamp	Displays the status of X-ray exposure. <ul style="list-style-type: none"> ▪ Green: Standby ▪ Yellow: In operation (X-ray on)
14	X-ray Generator	The X-ray tube where the X-ray is produced.
15	Vertical Frame	Holds the Rotating Unit. It can be controlled with the Column UP/DOWN switch.
16	Rotating Unit	Rotates around the patient's head while the image is being acquired. (Its movement is different according to the scan mode.)

1.6.1 Control Panel

No.	Buttons	Description
1	COLUMN UP/DOWN Button	Moves the Vertical Frame up or down. (For adjusting the height of the Chinrest)

1.6.2 Emergency Stop Switch

During operation, the following emergencies may occur:

- X-ray emission even after the **Exposure Switch** has been released
- Physical injury to the patient or damage to the equipment
- Other emergencies

If a problem occurs during image acquisition, press the red **Emergency Stop Switch** to immediately stop the moving parts and cut off all power to the equipment's electrical components. To restart the equipment, turn the **Emergency Stop Switch** clockwise until it pops up.

The **Emergency Stop Switch** is located under the bottom of the Handle Frame.

1.6.3 Exposure Switch

The **Exposure Switch** allows the operator to control image acquisition from outside of the X-ray room.

Press and hold the **Exposure Switch** down until acquisition is completed. Premature release of the **Exposure Switch** will abort image acquisition.

Pressing the **Exposure Switch** activates the LED indicator to turn yellow. This color indicates that the X-ray is being emitted.

IMPORTANT

The **Exposure Switch** is detachable. Ensure the **Exposure Switch** cable is not detached from the unit accidentally during the operation.

Keep vocal/visual contact with the patient during exposure. If any problem occurs during exposure, release the **Exposure Switch** immediately.

1.6.4 Chair purchase recommendations (The item that must be purchased separately)

This equipment is designed for the patient to sit in a chair and positioning. Therefore, Use this equipment after purchasing a separate chair. Here are some recommendations for choosing which chair to purchase.

Required chair specifications

- Chair without a backrest
- Up/Down available chair (requires height adjustment handle)
- When the chair is lowered to the lowest level or the upper part of the chair without vertical movement should be 450 mm or less from the ground
- When the chair is raised to the top level or the upper part of the chair without vertical movement should be 580 mm or more from the ground.
- Must have casters to facilitate chair movement

Optional chair specifications


- Chair with a backrest not too high
- Chair with an Anti-saddle-rotation function
- Chair with a footrest

NOTICE

- If a pediatric patient is an X-ray exposed while using a chair with a backrest if the shoulder line of the child is lower than the height of the backrest, use a child cushion to increase the sitting height of the pediatric patient (to prevent bump with the chair backrest and the equipment rotator part).

NOTICE

- For patients with reduced mobility or patients ((the elderly or children) who have difficulty fixing their body for a long time, if a chair with a 'backrest' is used, image quality degradation may occur due to the patient leaning on his/her back. Before using a chair with a backrest, it is required to check the appendix section.

This Page is Left Blank Intentionally


2. Getting Started

2.1 Turning on the Equipment

- Do not place the patient close to the equipment when it is being turned on. Doing so may cause physical injury to the patient and damage to the equipment.
- Do not operate the PC while the equipment is in operation. Doing so may cause an error in the equipment.

- Press the **Main Power Switch** that is located under the Handle frame to turn on the equipment.

NOTICE

Main Power Switch isolates its circuits electrically from the supply mains on all poles simultaneously.

- Make sure that the White LED light at the top of the equipment is on.

2.2 Running the Image Viewer (EzDent-i)

The Imaging Program is interfaced with **EzDent-i**, and the user can analyze the image acquired from the Console Software easily and rapidly. On your desktop, double-click the **EzDent-i** icon. The **EzDent-i** main window will be displayed.

NOTICE

For **vatech A9 (PHT-30CSS)** dental computed tomography X-ray system, both 3D viewer (**Ez3D-i**) and Console Software are being accessed through 2D viewer (**EzDent-i**) SW. 3D viewer and Console Software do not have an image storage capacity of their own, and both programs will not be able to keep patient information.

2.2.1 Creating a New Patient Record

To create a new patient record, follow the procedure outlined below:

1. Click the **PATIENT** tab and click the **Add Patient** icon from the main GUI window.

2. Enter the required patient information. **Chart Number**, **E-Mail address**, **First Name**, and **Last Name** are required fields that must be filled in. (The Chart Number is filled in automatically.)
3. Click the **Add** button to save the patient record.

2.2.2 Searching the Patient Records

You can search through the patient database using a patient's Chart Number, First Name, or Last Name.

1. Enter the Name or Chart Number of the patient to be searched on the **Patient Search** panel and then click the **Search** button. The information on the patient that fits the search condition appears.

Chart No.	Name	Date of Birth
20221031_102015	DiResta Jonathan	01/08/2022
20221031_102059	Wapelhorst Lloyd	01/01/2022
20221031_102144	Enter Manan	01/01/2022

NOTICE

Double-click the Keyboard icon to display the virtual keyboard. You may search for patient information using the virtual keyboard.

2. Double-click the patient information to see more details about the patient, as shown below.

2.3 Initiating the Console Software

NOTICE

For a new patient, the patient information needs to be registered first.

1. Search and select the patient to be captured.

Main Menu
EzDent-i

PATIENT ACQUISITION VIEWER CONSULT REPORT

SEARCH

RECENT PATIENTS LIST

Recently Acquired

Recently Viewed

DOCTOR

All

Chart No. 20221031_102059

Name Wapelhorst Lloyd

Gender/Age Male/0Y9M

Date of Birth 01/01/2022

More Details

Chart No.	Name	Date of Birth
20221031_102015	DiResta Jonathan	01/08/2022
20221031_102059	Wapelhorst Lloyd	01/01/2022

2. Click the **ACQUISITION** tab and select the imaging mode (CT, Panorama, or Cephalo).

Main Menu
EzDent-i

PATIENT ACQUISITION VIEWER CONSULT REPORT

SEARCH

CT

Panorama

Cephalo

Others

3. The Main Screen for the selected mode appears. From the Main Screen, you can configure the imaging parameter settings before acquiring an image.


NOTICE

Refer to the following **chapters (3~5)** for more information on image acquisition.

3. Acquiring PANO Images


3.1 Configuring Exposure Parameters

To acquire PANO Images, **2. Getting Started** must be completed first. If not, you must return to the **2. Getting Started** and finish the step first.

3. Acquiring PANO Images

1. Click the **PANO** button on the Main Screen.

NOTICE

The **CEPH** button exists only when the CEPH imaging program is included in the equipment.

2. Make an Arch Selection.

Arch Selection	Description
Narrow	Panoramic image of V-shaped palatal arches (Small number of adult females)
Normal	Panoramic image of normal adult palatal arches
Wide	Panoramic image of square-shaped palatal arches (some number of adult males)
Orthogonal	<p>The panoramic image where the x-ray angle enters vertically in between the teeth, so overlapping images is minimized.</p> <div> <div>NOTICE</div> <p>If Orthogonal Arch is selected, Bitewing examinations (Bitewing, Bitewing Incisor (Optional), Bitewing Right, Bitewing Left) are activated.</p> </div>
Child	Panoramic image of child palatal arches, more than 40% less X-ray dose than in Normal mode.

3. Select an Examination Program in the PANO Examination or Special Examination panel.

PANO Examination		
Right	Front	Left
Bitewing Right	Standard	Bitewing Left
Bitewing		

NOTICE

- To activate Bitewing examination options- **Bitewing**, **Bitewing Incisor** (Optional), **Bitewing Right**, **Bitewing Left**, select **Orthogonal** in the Arch Selection panel.

Arch Selection		
Narrow	Normal	Wide
Orthogonal	Child	

- When a Special Examination option is clicked, the **PANO Examination** panel is disabled. If you want to select a PANO Examination option, please conduct Arch selection again.

4. The Gender/Age group of the patient is selected automatically based on the patient information. If necessary, you can select the option manually.

NOTICE

Gender / Age Group		VATECH's Standard
Child		2 ~ 12 years of age
Adult	Man	> 12 years of age
	Woman	

- Select X-ray intensity.

NOTICE

Depending on the circumference of the patient's head, X-ray intensity may be classified as Hard, Normal, or Soft:
Soft ≤ Normal ≤ Hard

Age Group	Average Head Circumference (cm)	Range (cm)	X-ray Intensity
Child	53±3	>53±3	Hard
		53±3	Normal
		<53±3	Soft
Adult	56±3	>56±3	Hard
		56±3	Normal
		<56±3	Soft

- The values of tube voltage and current are configured automatically according to the patient's gender/age group and X-ray intensity. Click the **UP/DOWN** arrow to adjust kVp and mA. The dose is adjustable by ±1 kVp and ±0.1 mA, respectively.
- Click the **CONFIRM** button when the exposure parameter setting is completed.


CONFIRM

NOTICE

When you click **CONFIRM** button,

- The Rotating Unit will move to its initial scanning position.
- The Vertical Beam will be activated to make patient positioning easier.
- The DAP (Dose Area Product), Scan Time and Exposure Time will be displayed below the Scan Information window.

Scan Information

- Guide the patient to the equipment.

3.2 Patient Positioning

 WARNING	<ul style="list-style-type: none"> ▪ Have patience (especially pregnant women and children) wear a lead apron to protect themselves from residual radiation. ▪ Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.
IMPORTANT	<ul style="list-style-type: none"> ▪ Correct posture reduces the shadow cast by the patient's cervical spine and allows clear image acquisition. ▪ Metal implants or bridges may reduce the quality of the images. ▪ Be sure to adjust the laser beam correctly. Otherwise, the quality of images can be lower due to ghost images or expansion/reduction of the images.
IMPORTANT	<p>Since the patient needs to sit in a chair and position, the action that the user has to lower the upper body when positioning the patient is a necessary motion.</p>
NOTICE	<p>There is a horizontal beam laser for the guide that points the patient's shoulder in the rotator section. This laser allows you to proactively identify the possibility of moving the rotator and touching the patient's shoulder. If the laser is pointing at the patient's shoulder, it can be lowered further to prevent the rotator from touching the patient's shoulder.</p>
NOTICE	<p>The manufacturer recommends using a chair that can be adjusted to a height of at least 450 mm and up to 580 mm.</p>

Getting prepared

1. Let the patient remove all the metal objects (glasses, earrings, hairpins, braces, false teeth, etc.). Metal objects may induce ghost images and lower image quality.
2. Have the patient wear a lead apron to protect themselves from residual radiation.
3. Use the **Column UP/DOWN** button in the Handle frame or use the switch option to adjust the equipment to match the sitting height of the patient.

Chinrest position adjustment

Before patient alignment, adjust the Chinrest position to the PANO mode. (If the Chinrest position is already set to PANO mode, there is no need to adjust again.)

NOTICE

The installation angle of the Chinrest attachment part may be distorted by 2 degrees. This is intended for calibration purposes.


1. Push the Chinrest attachment part on the handle frame towards the machine until the end of the Chinrest attachment part is printed in "-PANO-" silk. Push until you hear a "click" sound.

3.2.1 PANO Examination Mode (Standard/Right/Left/Front/Orthogonal)

Normal Patient Positioning

1. Insert the **Normal Bite** and cover it with a **Sanitary Vinyl Cover**.

! CAUTION

- The Sanitary Vinyl Cover is for sole use only. It should be replaced for each patient. Be sure to use the approved vinyl cover.


- Clean the Chinrest and the Bite with non-alcohol based, non-corrosive cleaning solution and wipe with a dry towel before the next patient.

2. Guide the patient to a chair (sold separately) in the center of the equipment.
3. Guide the patient to sit in the center of the equipment


If a chair with a backrest is used, make sure that the patient's back is away from the backrest of the chair for proper positioning.

NOTICE


3. Acquiring PANO Images


4. Have the patient bite the Bite by mouth. (a little)

5. Position the patient's back and cervical spine as straight as possible.
6. Adjust the instrument to the patient's sitting height using the **Column UP/DOWN** on the Handle frame or the optional UP/DOWN switch while the patient holds the handle.


7. Close contact patient's chest to the handle frame.
8. Maintain the position described below.
 - Hold the handle firmly.
 - Make sure the patient's chest is in contact with the equipment.
 - Keep both shoulders parallel.
 - Straighten the patient's cervical spine and sit still.
 - Let the patient bite the Bite along its grooves with his/her front teeth.

9. Let the patient maintain the position described below.

- Close the mouth.
- Place the tongue to the roof of the mouth.
- Close the eyes.

Laser Beam Aligning


Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.

If the laser beam is not correctly positioned, there may be distortion, causing the image to be enlarged or reduced, or ghost shadows may occur and lower the image quality. Be sure to align the laser beam properly.

1. Align the Vertical Beam with the center of the face (Mid-sagittal Line). (It is to prevent the horizontal expansion of the image)
2. Align the Horizontal Beam in a straight line to the Frankfurt Line on the patient's face. Use the **Horizontal Beam** button on the control panel to position it. Make sure that the Horizontal Beam is aligned to the patient's face horizontally.

Vertical Beam /Mid-sagittal Line

Finishing Patient Positioning


1. Click the **READY** button. X-ray exposure has not started yet.

READY

2. Now go to **3.3 X-ray Exposure** to start the exposure.

Edentulous Patient Positioning

1. Insert the **Edentulous Bite**.


Clean the Edentulous Bite with non-alcohol based, non-corrosive cleaning solution and wipe with a dry towel before the next patient.

2. Guide the patient to a chair (sold separately) in the center of the equipment.
3. Guide the patient to sit in the center of the equipment

If a chair with a backrest is used, make sure that the patient's back is away from the backrest of the chair for proper positioning.

NOTICE

4. Position the patient's back and cervical spine as straight as possible.
5. Adjust the instrument to the patient's sitting height using the **Column UP/DOWN** on the Handle frame or the optional UP/DOWN switch while the patient holds the handle.
6. Close contact patient's chest to the handle frame.
7. Maintain the position described below.
 - Hold the handle firmly.
 - Make sure the patient's chest is in contact with the equipment.
 - Keep both shoulders parallel.
 - Straighten the patient's cervical spine and sit still.
 - Let the patient bite the Bite along its grooves with his/her front teeth.

8. Let the patient maintain the position described below.

Finishing Patient Positioning

1. Click the **READY** button. X-ray exposure has not started yet.

READY

2. Now go to **3.3 X-ray Exposure** to start the exposure.

3.2.2 SPECIAL Examination Mode (TMJ/Sinus)

[TMJ OPEN Mode (LAT)]

The TMJ Close image can be acquired after the TMJ Open image is acquired.


NOTICE

Steps for TMJ Mode

Patient positioning for TMJ Open → Laser Beam Aligning → X-ray Exposure → Patient positioning for TMJ Close → Laser Beam Aligning → X-ray Exposure

Patient Positioning

1. Remove the **Normal Chinrest** and insert the **Special Chinrest** into the equipment.
2. Insert the **Special Bite A** into the **Special Chinrest**.


! CAUTION

Clean the Chinrest and the Bite with non-alcohol based, non-corrosive cleaning solution and wipe with a dry towel before the next patient.

3. Guide the patient to a chair (sold separately) in the center of the equipment.
4. Guide the patient to sit in the center of the equipment

If a chair with a backrest is used, make sure that the patient's back is away from the backrest of the chair for proper positioning.

NOTICE

5. Position the patient's back and cervical spine as straight as possible.

3. Acquiring PANO Images


6. Adjust the instrument to the patient's sitting height using the **Column UP/DOWN** on the Handle frame or the optional UP/DOWN switch while the patient holds the handle.
7. Close contact patient's chest to the handle frame.
8. Guide the patient to press the base of the nose (acanthion point) against the Chinrest and tilt the head forward about 5°. At this point, make sure the patient's jaw does not touch the equipment.

IMPORTANT

- If the jaw touches the equipment, it is difficult to maintain the proper position to get good images.
- Be careful; the patient does not touch the equipment with his/her jaw.

9. Let the patient maintain the posture as follows:

- Open the mouth.
- Place the tongue to the roof of the mouth.
- Close the eyes.

NOTICE

- As shown in the picture, the support unit of the integrated Chinrest should touch the patient's acanthion point.
- Ask the patient to maintain his/her position until the operation is completed.

Laser Beam Aligning

WARNING


Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.

CAUTION

If the laser beam is not correctly positioned, there may be distortion, causing the image to be enlarged or reduced, or ghost shadows may occur and lower the image quality. Be sure to align the laser beam properly.

1. Align the Vertical Beam with the center of the face (Mid-sagittal Line). (It is to prevent the horizontal expansion of the image)
2. Align the Horizontal Beam in a straight line to the Frankfurt Line on the patient's face. Use the **Horizontal Beam** button on the control panel to position it. Make sure that the Horizontal Beam is aligned to the patient's face horizontally.

Finishing Patient Positioning

1. Click the **READY** button. X-ray exposure has not started yet.

READY

2. Now go to **3.3 X-ray** Exposure to start the exposure.

[TMJ CLOSE Mode (LAT) and SINUS (PA) Mode]

The TMJ Close image can be acquired after the TMJ Open image is acquired.

NOTICE

Steps for TMJ Mode


Patient positioning for TMJ Open → Laser Beam Aligning → X-ray Exposure → Patient positioning for TMJ Close → Laser Beam Aligning → X-ray Exposure

Patient Positioning

1. "Do you want to capture a TMJ Close image?" message will appear when the TMJ Open mode is completed. Press/Click **OK** button to begin TMJ Close mode.
2. Guide the patient to a chair (sold separately) in the center of the equipment.
3. Guide the patient to sit in the center of the equipment


If a chair with a backrest is used, make sure that the patient's back is away from the backrest of the chair for proper positioning.

NOTICE

4. Position the patient's back and cervical spine as straight as possible.
5. Adjust the instrument to the patient's sitting height using the **Column UP/DOWN** on the Handle frame or the optional UP/DOWN switch while the patient holds the handle.
6. Close contact patient's chest to the handle frame.
7. Guide the patient to sit in the center of the equipment and maintain the position described below.
 - Hold the handle firmly.
 - Make sure the patient's chest is in contact with the device.
 - Keep both shoulders parallel.
 - Straighten the patient's cervical spine and sit still.
 - Let the patient bite the Bite along its grooves with his/her front teeth.

8. Guide the patient to place the base of his/her nose (acanthion point) against the Chinrest and bend the head forward about 5°.
9. Let the patient maintain the posture as follows:
 - Close the mouth.
 - Place the tongue to the roof of the mouth.
 - Close the eyes.

NOTICE

- As shown in the picture, the support unit of the integrated Chinrest should touch the patient's acanthion point.
- Let the patient maintain his/her position until the operation is completed.

Laser Beam Aligning

This is the same as the one for TMJ Open mode.


Finishing Patient Positioning

This is the same as the one for TMJ Open mode.

[TMJ CLOSE Mode (PA) and Sinus Mode (LAT)]

Patient Positioning

1. Remove the **Normal Chinrest** and insert the **Special Chinrest** into the equipment.
2. Insert the **Special Bite B** into the **Special Chinrest**.


Clean the Chinrest and the Bite with non-alcohol based, non-corrosive cleaning solution and wipe with a dry towel before the next patient.

3. Guide the patient to a chair (sold separately) in the center of the equipment.
4. Guide the patient to sit in the center of the equipment

If a chair with a backrest is used, make sure that the patient's back is away from the backrest of the chair for proper positioning.

NOTICE

5. Position the patient's back and cervical spine as straight as possible.
6. Adjust the instrument to the patient's sitting height using the **Column UP/DOWN** on the Handle frame or the optional UP/DOWN switch while the patient holds the handle.
7. Close contact patient's chest to the handle frame.

8. Guide the patient to sit in the center of the equipment and maintain the position described below.
 - Hold the handle firmly.
 - Make sure the patient's chest is in contact with the device.
 - Keep both shoulders parallel.
 - Straighten the patient's cervical spine and sit still.
9. Guide the patient to press the base of the nose (acanthion point) against the Chinrest and tilt the head forward about 5°. At this point, make sure the patient's jaw does not touch the equipment.

IMPORTANT

- If the jaw touches the equipment, it is difficult to maintain the proper position to get good images.
- Be careful; the patient does not touch the equipment with his/her jaw.

10. Let the patient maintain the posture as follows:

- Close the mouth.
- Place the tongue to the roof of the mouth.
- Close the eyes.

NOTICE

- As shown in the picture, the support unit of the integrated Chinrest should touch the patient's acanthion point.
- Ask the patient to maintain his/her position until the operation is completed.

Laser Beam Aligning

Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.

If the laser beam is not correctly positioned, there may be distortion, causing the image to be enlarged or reduced, or ghost shadows may occur and lower the image quality. Be sure to align the Laser Beam properly.

1. Align the Vertical Beam with the center of the face (Mid-sagittal Line). (It is to prevent the horizontal expansion of the image)
2. Align the Horizontal Beam in a straight line to the Frankfurt Line on the patient's face. Use the **Horizontal Beam** button on the control panel to position it. Make sure that the Horizontal Beam is aligned to the patient's face horizontally.

Finishing Patient Positioning

1. Click the **READY** button. X-ray exposure has not started yet.

READY

2. Now go to **3.3 X-ray Exposure** to start the exposure.

3.3 X-ray Exposure

- If an emergency occurs during image acquisition, release the **Exposure Switch** to cease X-ray emission.
- The operator shall always observe the X-ray safety regulations applicable to his/her area during the operation of this equipment.

- The operator must always keep vocal/visual contact with the patient during the image acquisition process.
- Do not operate the PC during exposure. Doing so may cause the system to malfunction.

- Let the patient close the eyes during the operation.
- To acquire optimized images, instruct the patient to hold his/her breath and not to swallow. Also, do not let the patient move.

1. Get out of the X-ray room and close the door.

The operator must always keep vocal/visual contact with the patient during image acquisition.

2. Press and hold down the **Exposure Switch** until image acquisition is completed.

The image appears on the screen.

During X-ray exposure, the status appears as follows.

- The LED light of the **Exposure Switch** turns yellow.
- The LED light on the top of the equipment turns yellow.
- An alert sound comes out to indicate that X-ray emission is currently underway.
- On Console Software, the radiation mark turns Yellow.

3. Release the **Exposure Switch** when “Image capturing is completed” message appears on the screen.

3.4 Finishing the Scan

1. Guide the patient out of the equipment.
2. For **Normal Bite**, remove the **Sanitary Vinyl Cover** from the Bite.

3.5 Checking the Captured Images

Acquired images can be reconstructed and converted to DICOM format.

The exported images can be confirmed in **EzDent-i**.

NOTICE

Refer to the **EzDent-i User Manual** for more information.

1. The images are transferred to **EzDent-i** automatically.
2. The images are automatically saved if the automatic save option is configured as default. If it is not configured as default, click the **Save** button to save the images.
3. To check the image, double-click the one on the **Patient List**.


4. Acquiring CEPH Images (Optional)

4.1 Configuring Exposure Parameters

To acquire CEPH images, **2. Getting Started** must be completed first.


NOTICE

It is able to set the imaging parameters on the Console Software running on the PC. They are synchronized and display the same environmental settings.

4. Acquiring CEPH Images (Optional)


1. Click the **CEPH** button on the Main Screen.


NOTICE

The **CEPH** button exists only when the CEPH imaging program is included in the equipment.

2. Select an Image Option.

3. Select an examination program in the CEPH Examination panel.

4. The Gender/Age group of the patient is selected automatically based on the patient information. If necessary, you can select the option manually.

NOTICE

Gender / Age Group		VATECH's Standard
Child		2 ~ 12 years of age
Adult	Man	> 12 years of age
	Woman	

5. Select X-ray intensity.

NOTICE

Depending on the circumference of the patient's head, X-ray intensity may be classified as Hard, Normal, or Soft:
Soft ≤ Normal ≤ Hard

Age Group	Average Head Circumference (cm)	Range (cm)	X-ray Intensity
Child	53±3	>53±3	Hard
		53±3	Normal
		<53±3	Soft
Adult	56±3	>56±3	Hard
		56±3	Normal
		<56±3	Soft


6. The values of tube voltage and current are configured automatically according to the patient's gender/age group and X-ray intensity. Click the **UP/DOWN** arrow to adjust kVp and mA. The dose is adjustable by ±1 kVp and ±0.1 mA, respectively.
7. Click the **CONFIRM** button when the exposure parameter setting is completed.

CONFIRM

NOTICE

When you click **CONFIRM** button,

- The Rotating Unit will move to its initial scanning position.
- The Vertical Beam will be activated to make patient positioning easier.
- The **DAP** (Dose Area Product), Scan Time and Exposure Time will be displayed below the Scan Information window.

8. Guide the patient to the equipment.

4.2 Patient Positioning

 WARNING	<ul style="list-style-type: none"> ▪ Have patience (especially pregnant women and children) wear a lead apron to protect themselves from residual radiation. ▪ Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.
 CAUTION	Ensure that the Nasal Positioner is left unfolded, before adjusting the Ear Rods in the proper direction.
IMPORTANT	<ul style="list-style-type: none"> ▪ Correct posture reduces the shadow cast by the patient's cervical spine and allows clear image acquisition. ▪ Metal implants or bridges may reduce the quality of the images.
NOTICE	The manufacturer recommends using a chair that can be adjusted to a height of at least 450 mm and up to 580 mm.

Getting prepared

1. Let the patient remove all the metal objects (glasses, earrings, hairpins, braces, false teeth, etc.). Metal objects may induce ghost images and lower image quality.
2. Have the patient wear a lead apron to protect themselves from residual radiation.
3. Use the **Column UP/DOWN** button in the Handle frame or use the switch option to adjust the equipment to match the sitting height of the patient.

4.2.1 Lateral/Full Lateral (Optional) Mode

NOTICE


Correct posture reduces the shadow cast by the patient's cervical spine and allows clear image acquisition.

Patient Positioning

1. Turn the Nasal Positioner to the **Lateral** mode Positioning Marker, as shown below.

2. Leave enough space between the Ear Rods.


NOTICE

Use the Porion Position Reference Indicator that appears in the acquired image to confirm the location of Porion quickly.

3. Guide the patient to the CEPH unit.

4. Acquiring CEPH Images (Optional)

4. Instruct the patient to relax the neck and shoulders and sit upright.

5. Adjust the instrument to the patient's sitting height using the **Column UP/DOWN** on the Handle frame or the optional **UP/DOWN** switch.


Column UP/DOWN Button


WARNING

After adjusting the height of the column, align the Ear Rods and Nasal Positioner to the patient.

6. Align the Ear Rods into the patient's ears properly so that the head does not move during the operation. Moreover, aligning the Nasal Positioner with the patient's nasion by adjusting its height.

7. Align horizontally, so the patient's Frankfurt Line is parallel with the floor.

8. Direct the patient to swallow first before closing the mouth and to remain in his/her current position until image acquisition is completed.
9. Click the **READY** button. The x-ray exposure has not started yet.

READY

10. Now go to **4.3 X-ray Exposure** to start the exposure.


4.2.2 PA Mode

Patient Positioning


1. Turn the Nasal Positioner to the **PA/Waters' view/Carpus** mode Positioning Marker, as shown below.

2. Fold the Nasal Positioner up. The Nasal Positioner is not used in PA mode.
3. Guide the patient to the CEPH unit.

4. Ask the patient to sit upright towards the sensor. Make sure the patient's shoulders are flat and the neck is relaxed.


5. Adjust the instrument to the patient's sitting height using the **Column UP/DOWN** on the Handle frame or the optional **UP/DOWN** switch.

After adjusting the height of the column, align the Ear Rods to the patient.

6. During the operation, correctly align the Ear Rods to the patient's ears, so his/her head does not move.
7. Align horizontally, so the patient's Frankfurt Line is parallel with the floor.

8. Direct the patient to swallow first before closing his/her mouth and to remain in his/her current position until image acquisition is completed.
9. Click the **READY** button. The x-ray exposure has not started yet.

READY

10. Now go to **4.3 X-ray Exposure** to start the exposure.


4.2.3 SMV Mode

Patient Positioning

1. Turn the Nasal Positioner to the **SMV** mode Positioning Marker, as shown below.

2. Fold the Nasal Positioner up. The Nasal Positioner is not used in SMV mode.
3. Guide the patient to the CEPH unit.
4. Guide the patient to face the X-ray tube and sit upright.

5. Adjust the instrument to the patient's sitting height using the **Column UP/DOWN** on the Handle frame or the optional **UP/DOWN** switch.


WARNING

After adjusting the height of the column, align the Ear Rods to the patient.

6. During the operation, correctly align the Ear Rods to the patient's ears, so his/her head does not move.
7. Carefully tilt the patient's head back and adjust so his/her Frankfurt Line is vertical with the floor.

8. Direct the patient to swallow first before closing his/her mouth and to remain in his/her current position until image acquisition is completed.

9. Click the **READY** button. The x-ray exposure has not started yet.

READY

10. Now go to **4.3 X-ray Exposure** to start the exposure.

4.2.4 Waters' view Mode

Patient Positioning

1. Turn the Nasal Positioner to the **PA/Waters' view/Carpus** mode Positioning Marker, as shown below.


2. Fold the Nasal Positioner up. The Nasal Positioner is not used in Waters' view mode.
3. Guide the patient to the CEPH unit.
4. Ask the patient to sit upright facing the sensor. Make sure that the patient's shoulders are level and that his/her neck is relaxed.
5. Adjust the instrument to the patient's sitting height using the **Column UP/DOWN** on the Handle frame or the optional **UP/DOWN** switch.

After adjusting the height of the column, align the Ear Rods to the patient.

6. During the operation, correctly align the Ear Rods to the patient's ears, so his/her head does not move.

7. Direct the patient to swallow first before closing his/her mouth and guide the patient to bend the head backward 30° - 40° . Direct the patient to remain in the current position until image acquisition is completed.

8. Click the **READY** button. The x-ray exposure has not started yet.

READY

9. Now go to **4.3 X-ray Exposure** to start the exposure.

4.2.5 Carpus Mode


For Carpus Mode, install the Carpus Plate first before positioning the patient.

Installing the Carpus Plate

1. Turn the Nasal Positioner to the **PA/Waters' view/Carpus** mode Positioning Marker, as shown below.


2. Fold the Nasal Positioner up. The Nasal Positioner is not used in Carpus mode.
3. Fit the two ends of the Carpus Plate into the two holes of the CEPH unit as below.

4. Confirm that the Carpus Plate is safely mounted.

Patient Positioning

1. Let the patient put his/her right hand splayed on the Carpus Plate, as shown below. Make sure that the patient does not bend his/her fingers.

2. Ask the patient to close his/her eyes and sit still until the image acquisition is completed.
 3. Click the **READY** button. The x-ray exposure has not started yet.
- READY**
4. Now go to **4.3 X-ray Exposure** to start the exposure.

4.3 X-ray Exposure

- If an emergency occurs during image acquisition, release the **Exposure Switch** to cease X-ray emission.
- The operator shall always observe the X-ray safety regulations applicable to his/her area during the operation of this equipment.

- The operator must always keep vocal/visual contact with the patient during the image acquisition process.
- Do not operate the PC during exposure. Doing so may cause the system to malfunction.

- Let the patient close the eyes during the operation.
- To acquire optimized images, instruct the patient to hold his/her breath and not to swallow. Also, do not let the patient move.

1. Get out of the X-ray room and close the door.
2. Press and hold down the **Exposure Switch** until image acquisition is completed.

The image appears on the screen.

During X-ray exposure, the status appears as follows.

- The LED light of the **Exposure Switch** turns yellow.
- The LED light on the top of the equipment turns yellow.
- An alert sound comes out to indicate that X-ray emission is currently underway.
- On Console Software, the radiation mark turns Yellow.

3. Release the **Exposure Switch** when “Image capturing is completed” message appears on the screen.

4.4 Finishing the Scan

1. Leave enough space between the Ear Rods.
2. Fold the Nasal Positioner up in case it is unfolded.
3. Guide the patient out of the equipment.

4.5 Checking the Captured Images

Acquired images can be reconstructed and converted to DICOM format.

The exported images can be confirmed in **EzDent-i**.

NOTICE

Refer to the **EzDent-i User Manual** for more information.

1. The images are transferred to **EzDent-i** automatically.
2. The images are automatically saved if the automatic save option is configured as default. If it is not configured as default, click the **Save** button to save the images.
3. To check the image, double-click the one on the **Patient List**.

This Page is Left Blank Intentionally


5. Acquiring CBCT Images

5.1 Configuring Exposure Parameters


To acquire CBCT Images, **2. Getting Started** must be completed first.

NOTICE

You can set the imaging parameters on the Console Software running on the PC. They are synchronized and display the same environmental settings.

1. Click the **CBCT** button on the Main Screen.

NOTICE

The **CEPH** button exists only when the CEPH imaging program is included in the equipment.


NOTICE

Make sure FOV (Diameter x Height) is 8x8.

- ◆ Available options for each FOV are as below.

Available FOV (cm)	Vertical option	Horizontal option
8x8	Occlusion	Center


2. Select an Image Option.

3. Select a Voxel Size.

NOTICE

MAR (Metal Artifact Reduction) function is applied automatically if there are metal objects in the image. MAR may increase image reconstruction time.

4. The Gender/Age group of the patient is selected automatically based on the patient information. If necessary, you can select the option manually.

NOTICE

Gender / Age Group		VATECH's Standard
Child		2 ~ 12 years of age
Adult	Man	> 12 years of age
	Woman	

5. Select X-ray intensity.

NOTICE

Depending on the circumference of the patient's head, X-ray intensity may be classified as Hard, Normal, or Soft:
 $\text{Soft} \leq \text{Normal} \leq \text{Hard}$

Age Group	Average Head Circumference (cm)	Range (cm)	X-ray Intensity
Child	53±3	>53±3	Hard
		53±3	Normal
		<53±3	Soft
Adult	56±3	>56±3	Hard
		56±3	Normal
		<56±3	Soft

6. The values of tube voltage and current are configured automatically according to the patient's gender/age group and X-ray intensity. Click the **UP/DOWN** arrow to adjust kVp and mA. The dose is adjustable by ±1 kVp and ± 0.1 mA, respectively.

- Click the **CONFIRM** button when the exposure parameter setting is completed.

CONFIRM

NOTICE

When you click **CONFIRM** button,

- The Rotating Unit will move to its initial scanning position.
- The Vertical Beam will be activated to make patient positioning easier.
- The DAP (Dose Area Product), Scan Time and Exposure Time will be displayed below the Scan Information window.

Scan Information

- Guide the patient to the equipment.

5.2 Patient Positioning

WARNING

- Have patience (especially pregnant women and children) wear a lead apron to protect themselves from residual radiation.
- Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.

IMPORTANT

- Correct posture reduces the shadow cast by the patient's cervical spine and allows clear image acquisition.
- Metal implants or bridges may reduce the quality of the images.
- Be sure to adjust the laser beam correctly. Otherwise, the quality of images can be lower due to ghost images or expansion/reduction of the images.

NOTICE

There is a horizontal beam laser for the guide that points the patient's shoulder in the rotator section. This laser allows you to proactively identify the possibility of moving the rotator and touching the patient's shoulder. If the laser is pointing at the patient's shoulder, it can be lowered further to prevent the rotator from touching the patient's shoulder.

NOTICE

The manufacturer recommends using a chair that can be adjusted to a height of at least 450 mm and up to 580 mm.

IMPORTANT

Since the patient needs to sit in a chair and position, the action that the user has to lower the upper body when positioning the patient is a necessary motion.

Getting prepared

1. Let the patient remove all the metal objects (glasses, earrings, hairpins, braces, false teeth, etc.). Metal objects may induce ghost images and lower image quality.
2. Have the patient wear a lead apron to protect themselves from residual radiation.
3. Use the **Column UP/DOWN** button in the Handle frame or use the switch option to adjust the equipment to match the sitting height of the patient.

Chinrest position adjustment

Before patient alignment, adjust the Chinrest position to the CBCT mode. (If the Chinrest position is already set to CBCT mode, there is no need to adjust again.)

NOTICE

The installation angle of the Chinrest attachment part may be distorted by 2 degrees. This is intended for calibration purposes.


1. Pull the Chinrest attachment part on the handle frame from the machine until the end of the Chinrest attachment part is printed in "-CBCT-" silk. Pull until you hear a "click" sound.

Normal Patient Positioning

1. Insert the **Normal Bite** into the Normal Chinrest and cover it with a **Sanitary Vinyl Cover**.


- The Sanitary Vinyl Cover is for sole use only. It should be replaced for each patient. Be sure to use the approved vinyl cover.

- Clean the Chinrest and the Bite with non-alcohol based, non-corrosive cleaning solution and wipe with a dry towel before the next patient.


2. Guide the patient to a chair (sold separately) in the center of the equipment.
3. Guide the patient to sit in the center of the equipment

If a chair with a backrest is used, make sure that the patient's back is away from the backrest of the chair for proper positioning.


4. Position the patient's back and cervical spine as straight as possible.
5. Adjust the instrument to the patient's sitting height using the Column UP/DOWN on the Handle frame or the optional UP/DOWN switch while the patient holds the handle.

6. Close contact patient's chest to the handle frame.


7. Guide the patient to sit in the center of the equipment and maintain the position described below.

- Hold the handle firmly.
- Make sure the patient's chest is in contact with the device.
- Keep both shoulders parallel.
- Straighten the patient's cervical spine and sit still.
- Let the patient bite the Bite along its grooves with his/her front teeth.

8. Let the patient maintain the posture as follows:

- Close the mouth.
- Place the tongue to the roof of the mouth.
- Close the eyes.


Laser Beam Aligning – Normal

Be careful not to shine the Laser Beam directly into the person's eyes. Doing so may result in vision loss.

If the Laser Beam is not correctly positioned, there may be distortion, causing the image to be enlarged or reduced, or ghost shadows may occur and lower the image quality. Be sure to align Laser Beam properly.

Vertical Beam /Mid-sagittal Line

FOV 8x8 (cm) Center

NOTICE

This is a sample illustration for reference only. Actual FOV may vary from the image, as shown above.

1. Align the Vertical Beam with the center of the face (Mid-sagittal Line). (It is to prevent the horizontal expansion of the image)

Finishing Patient Positioning

1. After checking the positions of the patient and the Laser Beam, prevent the patient's head from moving.
2. Click the **READY** button. X-ray exposure has not started yet.

READY

- 9.
3. Now go to **9.5 X-ray Exposure** to start the exposure.

5.3 X-ray Exposure

- If an emergency occurs during image acquisition, release the **Exposure Switch** to cease X-ray emission.
- The operator shall always observe the X-ray safety regulations applicable to his/her area during the operation of this equipment.

- The operator must always keep vocal/visual contact with the patient during the image acquisition process.
- Do not operate the PC during exposure. Doing so may cause the system to malfunction.

- Let the patient close the eyes during the operation.
- To acquire optimized images, instruct the patient to hold his/her breath and not to swallow. Also, do not let the patient move.

1. Get out of the X-ray room and close the door.

The operator must always keep vocal/visual contact with the patient during image acquisition.

2. Press and hold down the **Exposure Switch** until image acquisition is completed.

The image appears on the screen.

During X-ray exposure, the status appears as follows.

- The LED light of the **Exposure Switch** turns yellow.
- The LED light on the top of the equipment turns yellow.
- An alert sound comes out to indicate that X-ray emission is currently underway.
- On Console Software, the radiation mark turns Yellow.

3. Release the **Exposure Switch** when “Image capturing is completed” message appears on the screen.

5.4 Finishing the Scan

1. Guide the patient out of the equipment.
2. For **Normal Bite**, remove the **Sanitary Vinyl Cover** from the Bite.

5.5 Checking the Captured Images

Acquired images can be reconstructed and converted to DICOM format.

The exported images can be confirmed in **EzDent-i**.

NOTICE

Refer to the **EzDent-i User Manual** for more information.






1. The images are transferred to **EzDent-i** automatically.
2. The images are automatically saved if the automatic save option is configured as default. If it is not configured as default, click the **Save** button to save the images.
3. To check the image, double-click the one on the **Patient List**.
Then, **Ez3D-i** automatically starts a 3D viewing of the saved image.

5.6 Image Comparison according to Patient Position – Incorrect Posture

5.6.1 Forward head posture (Turtle neck syndrome)

For patients with forward head posture, artifacts may occur as shown in the image below. (Cervical part Ghost image)





Position the patient's back and cervical spine in a position that is as straight as possible.

Normal C-shaped neck	Forward head posture (Turtle neck syndrome)
	
	
	 <p>Artifact Occurrence (Cervical Spine Part Ghost)</p>

5.6.2 Head-up posture

For patients with head-up posture, artifacts may occur as shown in the image below.
(arranged in a flat form image)





Make sure the patient does not up their head.

Normal C-shaped neck	Head-up posture
	
	 <p>Artifact Occurrence (arranged in a flat form image)</p>

5.6.3 Head-down posture

For patients with head-down posture, artifacts may occur as shown in the image below. (arranged in a V shape)




Make sure the patient does not down their head.

Normal C-shaped neck	Head-down posture
	
	 <p>Artifact occurrence (arranged in a V shape)</p>

5.6.4 Non-central location

For patients with non-central location posture, artifacts may occur as shown in the image below. (left/right asymmetry)

Guide the patient to be positioned in the central position.

Central location	Non-central location
	
	 <p>Artifact occurrence (left/right asymmetry)</p>

Copyright by © 2020 VATECH Co., Ltd.

All rights reserved.

The documentation, brand name, and logo used in this manual are copyrighted.

No part of this manual may be reproduced, transmitted, or transcribed without the expressed written permission of the manufacturer.

We reserve the right to make any alterations that may be required due to technical improvement. For the most current information, contact your VATECH representative.

Tel: (+82) 1588-9510

Email: gcs@vatech.co.kr

Website: www.vatech.com

Headquarters: 13, Samsung 1-ro 2-gil, Hwaseong-si, Gyeonggi-do, 18449, Korea

Factory: 13, Samsung 1-ro 2-gil, Hwaseong-si, Gyeonggi-do, 18449, Korea

The CE symbol grants this product compliance to the European Directive for Medical Devices 93/42/EEC as amended by 2007/47/EC as a class IIb device.

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Copyright de © 2020 VATECH Co., Ltd.