# EzRay Air™ Wall

Postal Code: 18449  
13, Samsung 1-ro 2-gil, Hwaseong-si, Gyeonggi-do, Korea  
[www.vatech.com](http://www.vatech.com)

**vatech**

EzRay Air™ Wall

User manual

For American region

# EzRay Air™ Wall

## User Manual

Model : VEX-S300W  
Version : 1.38

- English
- Français
- Português
- Español

For American regions

**vatech**

# ENGLISH

Notice.....	i
Conventions Used in this Manual .....	ii
1. System Overview .....	1
1.1 Indications for Use .....	1
1.2 Principles of Operation.....	1
1.3 Components.....	1
1.4 General View of the Equipment .....	3
2. Safety Instructions.....	9
2.1 General Safety Guidelines .....	9
2.2 Warnings and Safety Instructions .....	10
3. Operation.....	13
3.1 Power On/Off .....	13
3.2 Sleep Mode .....	14
3.3 Operation Mode .....	15
3.4 Positioning.....	20
3.5 Exposure .....	31
4. Troubleshooting.....	37
5. User Maintenance .....	39
6. Cleaning and Disinfection .....	41
6.1 Cleaning .....	42
6.2 Disinfection.....	42
Appendix .....	43
A.1 Using the Rotating Rectangular Cover .....	43

**This Page is Intentionally Left Blank**

## Notice

The **EzRay Air Wall (Model: VEX-S300W)** is an intra-oral dental X-ray system.

This manual contains descriptions, operational instructions, imaging procedures for the **EzRay Air Wall (Model: VEX-S300W)** dental X-ray system. It is recommended that you thoroughly familiarize yourself with this manual to make the most effective use of this equipment. Read and understand all cautions, safety messages and warnings in this manual.

Always keep this manual with the equipment and review the operation procedures and safety instructions if needed.

The illustrations/photos of the equipment in this manual are only for illustration purposes. Actual equipment may differ.

Due to continuous technological improvements, the manual may not contain the most updated information. For further information not covered in this manual, please contact us at:

**Phone:** (+82) 1588 9510

**E-mail:** gcs@vatech.co.kr

This document is originally written in English.

The **EzRay Air Wall (Model: VEX-S300W)** is referred to as **Equipment** or **System** in this manual.

**Manual Name:** EzRay Air Wall (Model: VEX-S300W) User Manual

**Version:** 1.38

**Publication Date:** 2023-06

## Conventions Used in this Manual

The following symbols are used throughout this manual. Make sure that you fully understand each symbol and follow the instructions accompanied.

To prevent physical injury and/or damage to the equipment, please observe all warnings and safety information included in this document.

 <b>WARNING</b>	<b>WARNING</b>	Indicates information that should be followed with the utmost care. Failure to comply with a warning may result in severe damage to the equipment or physical injury to the operator and/or patient.
 <b>CAUTION</b>	<b>CAUTION</b>	Indicates a situation that demands prompt and careful action, a specific remedy, or emergency attention.
<b>IMPORTANT</b>	<b>IMPORTANT</b>	Indicates a situation or action that could potentially cause problems to the equipment and/or its operation.
<b>NOTICE</b>	<b>NOTE</b>	Emphasizes important information or provides useful tips and hints.

# 1. System Overview

The **EzRay Air Wall (Model: VEX-S300W)** is an intra-oral dental X-ray system intended for intra-oral imaging. It consists of X-ray generator, X-ray controller, beam limiting device, operation panel and mechanical arm. The X-ray controller allows for accurate exposure control, and the adjustable mechanical arm allows for easy positioning. The system can be used with an imaging system.

## 1.1 Indications for Use

The **EzRay Air Wall (Model: VEX-S300W)** is an intra-oral dental X-ray system (extra-oral X-ray source system) intended for use by a trained and qualified dentist or dental technician for both adult and pediatric subjects for producing diagnostic dental radiographs for treatment of diseases of the teeth, jaw, and other oral structures using intra-oral image receptors.

## 1.2 Principles of Operation

X-rays are emitted when high voltage is supplied to the X-ray tube assembly which frees electrons from cathode. They hit anode to produce X-rays. The equipment acquires images by emitting X-rays continuously on human tooth.

## 1.3 Components


No.	Item	Standard	Option	Qty.
1	Power Box Assembly	•		1
2	Scissor Arm Assembly	•		1
3	X-ray Generator Assembly	•		1
4	User Manual	•		1
5	Installation Manual	•		1
6	Round Cover	•		1
7	Remote Exposure Switch*	•	•	1
8	Remote Exposure Switch (Doorbell type)	• (USA Only)	•	1
9	Remote Exposure Switch Cable (Doorbell type)	• (USA Only)	•	1
10	Door Interlock Cable	•	•	1
11	AC Power Cable		•	1
12	Rectangular Cover (2x3)		•	1
13	Rectangular Cover (4x3)	•	•	1

No.	Item	Standard	Option	Qty.
14	Double Stud	•	• (USA Only)	1
15	Horizontal Arm Assembly 450 mm		•	1
16	Horizontal Arm Assembly 600 mm		•	1
17	Horizontal Arm Assembly 900 mm		•	1
18	Scissor Arm Cover		•	1
19	Rotating rectangular cover (2x3)		•	1
20	Rotating rectangular cover (4x3)		•	1

\* The standard **Remote Exposure Switch** can be additionally provided as an optional item if needed. For details on the Remote Exposure Switch and Power Cable Connections, please see 'Section 3.5.1'.

## 1.4 General View of the Equipment

### Main Body

No.	Item	Description
1	Power Box Assembly	Controls the electrical power of the system.
2	Remote Exposure Switch and Power Cable Connectors	Connect the X-ray Exposure Switch cable and AC power cable (manufacturing option). Please refer to Section 3.5.1 Remote Exposure Switch and Power Cable Connections.
3	Main Power Switch	Main Power On/Off switch
4	Horizontal Arm	Connection arm between the power box assembly and the scissor arm assembly. (3 length options: 450 mm, 600 mm, 900 mm)
5	X-ray Beam Limiting Device	Limits the X-ray exposure area. Default type: Round Cone + Round Cover (FOV: Ø 6 cm)
6	Control Panel	Display for the X-ray exposure settings and operation conditions
7	X-ray Generator	Includes the X-ray tube and the high-voltage generator.
8	Handle	Grip the handle securely when using the system.
9	X-ray Generator Arm	Connection arm between the X-ray generator and the scissor arm assembly

## 1. System Overview

No.	Item	Description
10	Scissor Arm Assembly	Adjusts the length of the entire arm by folding and unfolding the scissor arm assembly. Connection arm assembly between the horizontal arm and the X-ray generator assembly.

## Control Panel

No.	Item	Description
1	Tube Voltage/Current Indicator	Indicates the tube voltage and tube current of the system.
2	Angle/Time Display	Displays the X-ray exposure time, error code, cooling time and exposure angle.
3	Adult/Child Selection	Indicates a patient type (adult or child).
4	Tooth Type Selection	Indicates a tooth type (incisor, canine, molar/premolar, bitewing).
5	X-ray Exposure Indicator	Indicates the X-ray exposure status. (Green: Ready / Yellow: X-ray ON)
6	SET Button (with LED lamp)	Resets the X-ray exposure angle. (SET button is pressed. → LED lamp flickers one time.)
7	Jog Dial	Turn the jog dial left (-) or right (+) to select X-ray exposure settings, press the jog dial to confirm the operating settings.
8	AUTO Button (with LED lamp)	Selects a tooth and exposure time automatically based on the X-ray exposure angle. (1. If the AUTO button is pressed, the LED lamp flickers one time. 2. If the AUTO Mode is ON, the LED lamp is ON.)
9	Speaker	Sound alarm for X-ray exposures

### Available Option Items

No.	Illustration/Photo	Option name	Usage	Material
1		Rectangular Cover 2x3 (FOV: 2x3 / 3x2 cm)	Used for limiting the X-ray exposure area by covering the X-ray Beam Limiting Device except for the 2x3 (3x2) rectangular area	ABS (Acrylonitrile butadiene styrene) copolymer
2		Rotating rectangular cover 4x3 (3x4) FOV: 4x3cm, 3x4cm (This adaptor can be used as both 4x3 and 3x4.)	Performs the same function as the Rectangular Cover. Moreover, Rectangular Cover rotates in 360 degrees and adjustable in 90-degree and 10-degree increments.	ABS (Acrylonitrile butadiene styrene) copolymer
3		Rotating rectangular cover 2x3 (3x2) FOV: 2x3cm, 3x2cm (This adaptor can be used as both 2x3 and 3x2.)	Performs the same function as the Rectangular Cover. Moreover, Rectangular Cover rotates in 360 degrees and adjustable in 90-degree and 10-degree increments.	ABS (Acrylonitrile butadiene styrene) copolymer
4		Remote Exposure Switch (if the Control Panel is installed to the wall framework)	Used to remotely perform exposure by pressing the Remote Exposure Switch	ABS (Acrylonitrile butadiene styrene) copolymer
5		Remote Exposure Switch (Doorbell type)	Used with the standard 'Remote Exposure Switch' (only if needed)	Steel (painted)
6		Remote Exposure Switch Cable (Doorbell type)	Used to connect the Remote Exposure Switch to the Remote Exposure Switch Connector on the bottom of the Power Box Assembly	PVC

No.	Illustration/Photo	Option name	Usage	Material
7		Door Interlock Cable	Used to connect the Door Interlock Switch to the Door Interlock Cable Connector on the bottom of the Power Box Assembly	PVC
8		Scissor Arm Cover	Used to cover the Scissor Arm	PVC
9		AC Power Cable	Used to connect the power (AC 220 V) to the Power Box Assembly	PVC
10		Double Stud (USA Only)	Used for fixing the power box to the wall (when installing the power box on the wooden wall in the North America regions)	ABS (Acrylonitrile butadiene styrene) copolymer + SPCC
11		Horizontal Arm 450 mm	Used as a connection arm (length: 450 mm) between the Power Box Assembly and the Scissor Arm Assembly	AL6061 + ABS (Acrylonitrile butadiene styrene)
12		Horizontal Arm 600 mm	Used as a connection arm (length: 600 mm) between the Power Box Assembly and the Scissor Arm Assembly	AL6061 + ABS (Acrylonitrile butadiene styrene)
13		Horizontal Arm 900 mm	Used as a connection arm (length: 900 mm) between the Power Box Assembly and the Scissor Arm Assembly	AL6061 + ABS (Acrylonitrile butadiene styrene)

**NOTICE**

Refer to the **Appendix A.1 'Using the Rotating Rectangular Cover'** for instructions for use.

**This page intentionally left blank**

## 2. Safety Instructions

### 2.1 General Safety Guidelines

- Mode of operation: Continuous operation with cyclic loading—This equipment needs a rest time of at least 60 times the exposure time before starting the next exposure.
- This equipment is designed and manufactured to ensure maximum safety of operation. Operate and maintain it in strict compliance with the safety precautions and operating instructions contained in this manual.
- This equipment must only be operated by legally qualified persons.
- Observe all local fire regulations. Always keep a fire extinguisher near the equipment.
- The equipment must be installed, maintained, and serviced by qualified service personnel according to the procedures and preventive maintenance schedules.
- Ensure that the on/off switch is set to off when the equipment is not in use.
- Always disconnect the power supply before cleaning the equipment.
- Don't keep the equipment or its parts in a humid place or near a liquid substance.
- Avoid placing the equipment near chemical storage and gas-filled storage facilities.

## 2.2 Warnings and Safety Instructions

 <b>WARNING</b>	<p>This X-ray unit may be dangerous to patient and operator unless safe exposure factors, operating instructions and maintenance schedules are observed.</p> <p>It is important to read this user manual carefully and strictly abide by all warnings and cautions stated within it.</p>
 <b>WARNING</b>	<p>To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.</p>
 <b>WARNING</b>	<p>Since rules and regulations concerning radiation safety differ between countries, it is the responsibility of the owner and/or operator of this equipment to comply with all applicable rules and regulations concerning radiation safety and protection in their area.</p>

- DO NOT open or remove the cover panels on this equipment.
- Never expose this equipment to liquids, mists or sprays. Exposing this equipment to liquids may cause an electrical shock or otherwise damage the system.
- DO NOT use spray cleaners on this equipment, as this could cause a fire.
- Never use this equipment in an environment that is susceptible to explosion.
- DO NOT place flammable materials near this equipment.
- Never touch the patient while also touching the SIP/SOP connectors.
- Never try to modify this equipment, including the wires or cables. Modifying this equipment may damage it beyond repair.
- We recommend that the patient and the operator wear protective lead-lined aprons, unless other Radiation Protection Protocols apply locally.
- Children and pregnant women must consult with a doctor before X-ray exposure.
- Serious dangers may occur from electromagnetic interference (i.e. noise) between other equipment in the area during specific examinations or medical treatment.

**Radiation Safety**

When using the equipment, it is recommended that all users comply with the following radiation safety guidance for the safety of the users and the patients.


- This equipment should be operated by a trained and qualified dentist or a dental technician in a controlled environment.
- All users and patients should wear protective equipment, such as a lead apron, thyroid collar, etc.
- Pregnant women should not be exposed to X-rays unless it is strictly necessary.
- All users should comply with the Radiation Protection Policies established by the government.
- Any person or organization who installs an external Door Interlock Switch is responsible for ensuring that the Switch has a radiation indicator or an equivalent alarm system to show the state of current.
- This equipment should be operated at least 2 m away from the operator when operated in or outside of the room.

**This page intentionally left blank**

### 3. Operation

#### 3.1 Power On/Off

1. Turn on the system referring to the following figure and table.

Main Power On/Off switch

2. The following displays light up. For further information, see 'Control Panel' on *page 5*.

- 1) Current Angle/Time Display
- 2) Adult/Child Selection Display
- 3) Tooth Type Selection Display
- 4) X-ray Exposure Indicator

### 3.2 Sleep Mode

If the SET button or AUTO button is not pressed or the jog dial is not turned for 5 minutes, the DISPLAY turns OFF.

In the Sleep mode, the LEDs of the SET button and AUTO button blink at 1 second intervals.

#### NOTICE

To off the Sleep mode,

- press the SET button or AUTO button.
- press or turn the jog dial.
- move header up or down

### 3.3 Operation Mode

This system can be operated with Manual Mode and Auto Mode, and you can set up the mode by using the **AUTO** button.

#### Manual Mode

1. To start the Manual Mode, check if the lamp under the **AUTO** button is turned off.

### 3. Operation

- When the tooth type selection area flickers, turn the jog dial to select the tooth type. To see the Control Panel before and after selection, refer to the figures below.

**Before tooth type selection**

**After tooth type selection**


**Tooth Type**

Symbol	Type
	Incisor
	Canine
	Molar/Premolar
	Bitewing

3. After tooth type selection, patient type should be selected. When the Adult/Child selection area flickers, turn the jog dial to select the patient type. To see the Control Panel after selection, refer to the figure below.

After patient type selection

Patient Type

Symbol	Type
	Adult
	Child

**NOTICE**

After the tooth type and patient type are selected, exposure time is automatically displayed.

4. If you want to change the exposure time, turn the jog dial to adjust the exposure time from 0.05 to 0.5 s. (increments: 0.01 s)

**NOTICE**

If you press the jog dial after adjusting the exposure time in Manual Mode, the exposure time is returned to the default setting.

To save the exposure time as default in Manual Mode, press and hold the jog dial for about 3 seconds.

#### **Auto Mode**

1. When the Auto Mode is turned on by pressing **AUTO**, the default angle is displayed as shown in the following figure.

#### **NOTICE**

To set the starting point during exposure, press **SET**.

2. Position the system to the teeth to perform exposure. To see how to perform the patient positioning, refer to Section 3.4 Positioning. (The following figure is an example of the maxilla incisor.)

3. When tooth types are selected, the exposure angles are automatically set according to the tooth type. To check the default exposure angles, refer to the following table.

Tooth Type	Angle of Inclination
Incisor	Maxilla: +40° ~ +50°
	Mandible: -22° ~ -28°
Canine	Maxilla: +40° ~ +50°
	Mandible: -17° ~ -23°
Molar/Premolar	Maxilla: +25° ~ +35°
	Mandible: -2° ~ -8°
Bitewing	+3° ~ +12°

### NOTICE

Since the angles of inclination for the maxilla incisor and canine are same, the exposure time of the canine is applied for both the maxilla incisor and canine.

Refer to the following figure to see the angle for molar/premolar.


### 3.4 Positioning

#### Positioning the Patient

To obtain high-quality intra-oral radiography with maximum details, take extra care in all steps of the radiography process: positioning the patient and the X-ray imaging system; exposing the intra-oral sensor.


1. Place a protective lead apron on the patient's chest.
2. Have the patient sit on the chair with the sagittal plane vertical.
  - For radiography of the upper maxillary, the Frankfort Plane must be horizontal.

- For radiography of the lower maxillary, the Occlusal Plane must be horizontal.


3. Place the tube head cone on the area you want to take an image.


For further information about the patient positioning and beam angle for each mode, refer to the following **Positioning Instructions**.

### Positioning Instructions

**Paralleling technique:** The sensor is placed in a holder which is used to align the sensor parallel to the long axis of the teeth.

**Bisected angle technique:** The patient holds the sensor in place with his/her finger. The X-ray beam is directed perpendicularly towards an imaginary line, which bisects the angle between the sensor plane and the long axis of the tooth.

Position the tube head to the patient using the accepted standard positioning procedures.


Here are the specific angulations and directions for the tube head in order to take the best images of a particular tooth (i.e. **Bisected angle technique**).

Position the receptor carefully not to damage the soft tissue of the patient's intra-oral area.

- **Maxillary Incisor**


X-ray beam is directed downward at 45°.

Teeth		Angle of inclination
Incisor	Maxilla	+45°

- **Mandibular Incisor**


X-ray beam is directed upward at 25°.

Teeth		Angle of inclination
Incisor	Mandible	-25°

- **Maxillary Canine**


X-ray beam is directed downward at 45°.

Teeth		Angle of inclination
Canine	Maxilla	+45°

- **Mandibular Canine**

X-ray beam is directed upward at 20°.

Teeth		Angle of inclination
Canine	Mandible	-20°

- **Maxillary Molar and Premolar**


X-ray beam is directed downward at 30°.

Teeth		Angle of inclination
Molar and Premolar	Maxilla	+30°

- **Mandibular Molar and Premolar**

X-ray beam is directed upward at 5°.

Teeth		Angle of inclination
Molar and Premolar	Mandible	-5°

- **Bitewing**

For a bitewing exposure, the patient closes their teeth during exposure on the sensor holder.

X-ray beam is directed downward at  $5^{\circ} \sim 8^{\circ}$ .

Teeth	Angle of inclination
Bitewing exposure	$+5^{\circ} \sim +8^{\circ}$

### Positioning the Imaging Sensor


To ensure image quality, the digital imaging sensor must be positioned properly (for information about the proper placement of the imaging sensor, please refer to 'Positioning Instructions').

- Failure to position the imaging sensor properly can result in errors on the radiograph, such as distorted teeth and roots, elongation, magnification, and overlapping contacts.


#### **NOTICE**

The paralleling technique generally reduces the risk of such errors, but if you position the sensor improperly, angulation errors may occur (angulation of the sensor to the tooth itself).

- Failure to align the imaging sensor with the exit pattern of the X-ray beam can result in cone-cuts on the radiograph. The cone-cuts are clear areas that are shown on the radiograph when part of the radiograph is not exposed to radiation. Please refer to the following figure as an example of cone-cuts.

The following figure indicates how the cone-cut occurred by showing the position of the imaging sensor and the radiation area.

To ensure proper alignment between the imaging sensor and the X-ray beam, it is recommended to use a PID (Position Indicating Device).

When using the PID, the exit pattern of the X-ray device should be aligned perpendicular to the target receptor.

---

**NOTICE**

Once the PID is properly aligned, instruct the patient not to move.


---

## 3.5 Exposure

**IMPORTANT**

The operator **MUST** instruct the patient to refrain from moving during the entire exposure.

1. Instruct the patient not to move.
2. Press and hold the **Remote Exposure Switch** for exposure duration. The **Remote Exposure Switch** allows the operator to control image acquisition from outside of the X-ray room. Pressing the **Remote Exposure Switch** activates the X-ray Exposure Indicator to turn yellow. This color indicates that the X-ray is being emitted.

**IMPORTANT**

Press and hold the **Remote Exposure Switch** as long as the acoustic signal can be heard. Otherwise the exposure will be faulty and there will be an error message on the Control Panel.


3. While X-ray is being exposed, the X-ray Exposure Indicators on the **Control Panel** and the **Remote Exposure Switch** are turned on and an audible sound is produced. Keep pressing until the X-ray Exposure Indicator lights are turned off and the audible sound stops.

<b>IMPORTANT</b>	The <b>Remote Exposure Switch</b> is detachable. Ensure that the <b>Remote Exposure Switch</b> cable is not detached out from the unit accidentally during the operation.
<b>IMPORTANT</b>	Keep vocal/visual contact with the patient during exposure. If any problem occurs during exposure, release the <b>Remote Exposure Switch</b> immediately.
<b>IMPORTANT</b>	If the <b>Door Interlock Switch (option)</b> has been installed, X-ray exposure will be stopped once the door of the X-ray room is opened.
<b>NOTICE</b>	As described in step 2 and 3 above, the X-ray Exposure Indicator is included both on the <b>Control Panel</b> and the <b>Remote Exposure Switch</b> .

### 3.5.1 Remote Exposure Switch and Power Cable Connections

There are four connectors on the bottom of the Power Box Assembly as described below.

- **Connector 1, 2, 3:** Remote Exposure Switch Connectors
- **Connector 4:** AC Power Cable Connector


The standard and optional **Remote Exposure Switches**, and the **Door Interlock Cable** are connected to the **Connector 1, 2, and 3**. Please check all seven possible options in the table below.

#### IMPORTANT

Before connecting the cables, ensure that the DIP switch setting is set according to the instructions in the following table.

Refer to the location of the DIP switch in the following figure.

Option No.	Description	Connector 1	Connector 2	Connector 3	Connector 4
Option 1	Press (1).		N/A	N/A	AC Power Cable (Option)
		(1) Remote Exposure Switch (Standard)			
Option 2	Press (1) and (2) at the same time.		N/A		AC Power Cable (Option)
		(1) Remote Exposure Switch (Standard)		(2) Remote Exposure Switch (Option)	
Option 3	Press (1) and (2) at the same time.		N/A		AC Power Cable (Option)
		(1) Remote Exposure Switch (Standard)		(2) Remote Exposure Switch (Doorbell type) (Option)	
Option 4	Press (1).		N/A		AC Power Cable (Option)
		(1) Remote Exposure Switch (Standard)		(2) Door Interlock Cable (Option)	
Option 5	Press (1) or (2).			N/A	AC Power Cable (Option)
		(1) Remote Exposure Switch (Standard)	(2) Remote Exposure Switch (Option)		

Option No.	Description	Connector 1	Connector 2	Connector 3	Connector 4
Option 6	Press (1) <b>or</b> (2).			N/A	AC Power Cable <b>(Option)</b>
		(1) Remote Exposure Switch <b>(Standard)</b>	(2) Remote Exposure Switch (Doorbell type) <b>(Option)</b>		
Option 7	Press (2) <b>and</b> (3) at the same time. (1): Not used in this option.				AC Power Cable <b>(Option)</b>
		(1) Remote Exposure Switch <b>(Standard)</b>	(3) Remote Exposure Switch (Doorbell type) <b>(Option)</b>	(2) Remote Exposure Switch (Doorbell type) <b>(Option)</b>	

**This page intentionally left blank**

## 4. Troubleshooting

In instances of abnormal operation, alarm/error messages will be displayed on the Control Panel. If a problem persists, please request assistance from the customer support information services.

### Alarm/Error Messages

#### NOTICE

**A.0X:** A problem occurred, and the system performs the correction automatically. This alarm clears after the correction is completed.

**E.0X:** An error occurred. Turn the power off, and then turn it back on. If the error persists, contact your Service Representative.

Error Code	Check Parameter	Description
E.02	X-ray Generator	Error related to X-ray exposure is not possible to exposure X-ray in the state while "E.02", "E.03", "E.04", "E.05" where power is maintained.  After X-ray exposure related error occurs, when the equipment is turned off and turned on, X-ray exposure is normally performed.
E.03		
E.04		
E.05		
A.06	X-ray Generator	The system needs cooling time due to continuous operation. (Duty Cycle: <b>1:60</b> ) This alarm clears when the system temperature goes down to normal.
A.07	System	The Remote Exposure Switch has been held down. Release the switch.
A.08	User	The exposure button has been pressed and released before the X-ray exposure finished with the default exposure time. Press and hold the exposure button for the duration of the exposure time.

**Troubleshooting**

Problem	Cause	Solution
Equipment is not turned on.	Power button is not turned on properly.	Turn the equipment power switch off and turn it back on.
Control Panel is not turned on.	Defective main board	Contact your Service Representative.
	Internal cable disconnected	Contact your Service Representative.
No X-ray emission	Generator is cooling.	Wait for the cooling time. (Duty Cycle: <b>1:60</b> )
	Remote Exposure Switch is pressed, but X-ray exposure does not work.	Contact your Service Representative.
	Internal cable disconnected	Contact your Service Representative.
	Defective generator	Contact your Service Representative.
	Tube lifecycle termination	Contact your Service Representative.
X-ray emission works, but exposure is too light or completely white.	Equipment has been positioned incorrectly.	Adjust the position of the equipment.
	Exposure time is too short.	Increase the exposure time.
	Receptor is facing the wrong way.	Reposition the receptor.
X-ray emission works, but exposure is too dark.	Exposure time is too long.	Decrease the exposure time.

## 5. User Maintenance

To ensure that the equipment operates at maximum efficiency, it is recommended to perform the following procedures for the customer's routine maintenance activities.

### Maintenance Task Checklist

Always turn off the equipment before performing any maintenance.

Tasks	Period
Sterilize all components that come into contact with the patient and operator by using an alcohol-based solution.	Daily
Wipe the outer covers of the equipment with a dry cloth at the end of each day's operation.	Daily
Ensure that the main power button has been turned off after using the equipment.	Daily
Ensure the audible signal is audible and the X-ray exposure light is visible when you perform an exposure.	Daily
Ensure that the wall framework is securely attached to the wall.	Daily
Ensure that the yellow (exposure) indicator light turns on when the exposure button is pressed.	Daily
Ensure that the power cable does not have cuts or abrasions.	Monthly
Ensure that all cable connectors are not mechanically defective.	Monthly
Ensure that the control panel has no defects.	Monthly
Ensure that the power cable does not have cuts or abrasions.	Monthly
Ensure that all visible labels are intact and legible.	Monthly

DO NOT use detergents or solvents to clean the outer covers of the equipment.

If any defects are found, do not operate the equipment since it has to be handled by a qualified person. Contact your Service Representative.

**This page intentionally left blank**

## 6. Cleaning and Disinfection

The liquid can cause damage to the equipment. When cleaning or disinfecting, liquids may enter the equipment or the release button via the ventilation slots.

- Do not spray the equipment with Cleaners or disinfectant agents. Apply the Cleaners or disinfectant agents to a clean cloth and wipe it.
- Make sure that no liquids run along the surface into the ventilation slots or release button.
- Remove any soiling with a soft, wet, lint-free cloth.
- Please follow the Cleaners or disinfectant agents' instructions for use.
- When cleaning or disinfect the surfaces, always disconnect the equipment from mains.
- Do not use spray cleaner or disinfectant directly onto the equipment, as this could cause a fire.

- Cleaners or disinfectant agents may contain powerful ingredients. Unsuitable cleaning and disinfectant agents are detrimental to health and attack the surface of the equipment.
- Do not use cleaners or disinfectant agents containing Phenol, acetic acid, peroxide, or other oxygen splitting agents, sodium hypochlorite, isopropyl alcohol (2-propanol, isopropanol) or iodine- splitting agents.
- Comply with the specifications contained in the operating instructions of the Cleaners or disinfectant agents.
- Wear safety gloves.

### 6.1 Cleaning

- The equipment surfaces can be cleaned with a soft cloth dampened in a non-alcohol based, non-corrosive cleaning solution. Do not use sponges or, in any case, any material that can be reused.

### 6.2 Disinfection

- Use only disinfectants that comply with the valid requirements of the respective national regulatory body or whose bactericidal, fungicidal, and virucidal properties have been verifiably tested and approved accordingly.
- Sterilization and disinfection should be performed thoroughly for items that have been in frequent contact with patients and operators.
- Do not use UV systems to disinfect the equipment, as exposed parts of the device can turn yellow or discolor.
- The use of unsuitable Cleaners or disinfectant agents and methods can damage the equipment and accessories. Only use the Cleaners or disinfectant agents specified or approved by VATECH
- The following Cleaners or disinfectant agents have been evaluated for safe use on the surfaces.
- Never combine products or with liquids other than the products listed above.
- Damages to surfaces and materials due to the use of different products cannot be excluded even if they are not included in the exceptions mentioned above.
- Use a non-alcoholic chlorine dioxide-based disinfectant.

#### Sample


## Appendix

### A.1 Using the Rotating Rectangular Cover


Rotating rectangular cover rotates in 360 degrees.

Rotating rectangular cover consists of the Outer Ring and the Inner Ring.


The Outer Ring rotates in 90-degree increments.

**IMPORTANT**

The Inner Ring rotates in 10-degree increments.

1. Assemble the rotating rectangular cover to the Cone of the Main body.

**NOTICE**

For assembling the rotating rectangular cover to the Cone, turn the rotating rectangular cover 45 degrees and check the Click sound to make sure the assembly is complete.

2. Turn the Outer Ring of the rotating rectangular cover to adjust the angle roughly.
3. Turn the Inner Ring of the rotating rectangular cover to fine-tune the angle.

**NOTICE**

Disassembly method of the rotating rectangular cover is the same as the assembly method. Refer to the NOTE mentioned below stage 1.

**Copyright by © 2017 VATECH Co., Ltd.**

All rights reserved.

The documentation, brand name and logo used in this manual are copyrighted.

No part of this manual may be reproduced, transmitted, or transcribed without the expressed written permission of the manufacturer.

We reserve the right to make any alterations which may be required due to technical improvements. For the most current information, contact your VATECH representative.

**Manufactured by VATECH Co., Ltd.**

**Tel:** (+82) 1588 9510

**Email:** gcs@vatech.co.kr

**Website:** www.vatech.com

**HQ Address:** 13, Samsung 1-ro 2-gil, Hwaseong-si, Gyeonggi-do, 18449, Korea

**Factory Address:** 13, Samsung 1-ro 2-gil, Hwaseong-si, Gyeonggi-do, 18449, Korea

The CE symbol grants this product compliance to the European Directive for Medical Devices 93/42/EEC as amended by 2007/47/EC as a class IIb device.

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**Copyright © 2017 VATECH Co., Ltd.**
**Fabricado por VATECH Co., Ltd.**